

# CLEANPatch-V<sup>®</sup>

## PRODUCT INFORMATION

Product Code	Clean Patch-V Size	Application
VP-S1	<u>Small Kit:</u>	<u>Damage up to:</u>
	2 x Strips 44mm x 100mm 1 x Rectangle 145mm x 200mm	5cm 15cm
VP-L1	<u>Large Kit:</u>	<u>Damage up to:</u>
	7 x Strips 44mm x 100mm 3 x Rectangle 145mm x 200mm	5cm 15cm

# CLEANPatch-V<sup>®</sup>

*Medical Surface Repair patch for vinyl upholstery.  
Restores damaged surfaces to an intact and hygienic state.*



For more information, please visit [www.cleanpatch.ca](http://www.cleanpatch.ca)

Manufactured by  
Surface Medical Inc.  
Calgary, AB, Canada T2S 1W2  
1-888-623-7085  
[info@surfacemedical.ca](mailto:info@surfacemedical.ca)

© Surface Medical Inc. 2018

Distributed in SA by



**SURGICAL AND MEDICAL SUPPLIES PTY LTD**  
16 Kensington Road Rose Park SA 5067  
Telephone: (08) 8332 1666 Fax: (08) 8332 5132  
Email: [sales@surmed.com.au](mailto:sales@surmed.com.au)  
Website: [www.surmed.com.au](http://www.surmed.com.au)

Contact  
Chris Lee, Product Specialist  
0476 777 767

## THE PROBLEM

Vinyl surfaces are prevalent in healthcare facilities, from wheelchairs to examination tables to geriatric chairs. These surfaces are frequently damaged and may harbor pathogens that contribute to Healthcare Associated Infections.

Until now, there was no simple solution to enhance patient and resident safety while reducing costs.

## THE SOLUTION IS HERE . . .

## PREVIOUS ALTERNATIVES

Replace or Recover	✓ Validated ✗ Unaffordable ✗ Inconvenient
Tape the damage	✓ Inexpensive ✗ Not Validated ✗ Not Durable
Do Nothing	✗ Risk Patient Safety



## WHY CLEANPATCH®-V?

CleanPatch-V is the first medical surface repair patch for vinyl surfaces, including wheelchairs, examination tables and geriatric chairs. It restores damaged surfaces to an intact and hygienic state.

**Clinically Tested**

- Durable and impervious to fluids under rigorous conditions
- Fully cleanable with a wide variety of hospital disinfectants
- Validated in healthcare settings

**Safe and Simple**

- Early intervention reduces the risk of cross contamination
- Made of healthcare grade materials
- Simple peel and stick application, adheres in seconds

**Cost Effective**

- Extends the life of a variety of vinyl surfaces
- Proactively improves patient safety
- Significant savings from vinyl repair versus recover or replacement of surfaces

## CleanPatch®-V should be used appropriately as part of a regular equipment maintenance program

	✓ Use CleanPatch-V	✗ Do Not Use CleanPatch-V
<b>Equipment Condition</b>	Usable within expected lifespan No signs of fluid ingress	No longer suitable for use Surface is stained or warped
<b>Damage Size and Location</b>	Less than 6 inches (15 cm) Flat areas or over edges (2-D) Not over seams	More than 6 inches (15 cm) Sharp corners (3-D) Over seam
<b>Application</b>	One patch per damaged area	Layering multiple patches is not recommended
<b>Examples</b>	 	 

To ensure product integrity, do not cut CleanPatch



### References:

FDA Safety Communication, April 2013

Oliviera, AC. et al. Journal of Nursing UFPE on line, Jan 2013

Otter JA, et al. Infection Control and Epidemiology, Vol 32 (7) July 2011

Viana, REI H. et al. American Journal of Infection Control, 1;44(4) April 2015