

SIGVARIS

NEW SOFT SILHOUETTE LEGGINGS COMPRESSION WITH BENEFITS

SOFT SILHOUETTE LEGGINGS feature 15–20mmHg of compression at the ankle. The level of compression decreases up the leg to improve circulation, revitalize legs and help prevent varicose and spider veins.

FEATURES

Stylish subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics and shirts

Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort

Advanced soft stretch microfiber fabric shapes legs and smooths skin

BEST SUITED FOR

Tired, heavy, restless, aching legs

Travelers

Occupational leg health risks

FIBERS

83% Nylon

17% Spandex

Latex-Free



BLACK
(99)

To order:
sales@surmed.com.au | (08) 8332 1666

SIGVARIS

NEW SOFT SILHOUETTE LEGGINGS COMPRESSION WITH BENEFITS



15–20mmHg Legging
Size B in Black
170LB99




Non-binding ankle cuff stays
comfortably in place for all day wear

The perfect pair of
leggings to keep your
legs looking great and
feeling great every day!

Product Code Example: For product at far left: 170LB99	Series	Product Style	Size	Color
	170	L	B	99

SIZING CHART



GUIDELINES FOR ANKLE CIRCUMFERENCE

SIZE A: 7"–9", SIZE B & C: 8"–10", SIZE D: 8.5"–11.5"