

SIGVARIS



FALL INTO
CLASSIC STYLE
ARGYLE
THE NEWEST MICROFIBER SHADES





MICROFIBER SHADES

GRADUATED COMPRESSION SOCKS WITH PERSONALITY

SIGVARIS MICROFIBER SHADES stripes and heathers have been best-sellers since their introduction, and now there is a new argyle pattern in four colors to entice your customers. Also, all of our MICROFIBER SHADES are now available in 20–30mmHg (six measured sizes) and 15–20mmHg (three easy-to-shop sizes). With our MICROFIBER SHADES line, SIGVARIS continues to prove that compression can be fashionable and functional.

FEATURES

- 100% microfiber next to skin for softness and excellent moisture management
- Double-covered inlay yarns for easy donning and doffing

BEST SUITED FOR

- Management of venous insufficiency
- Long-distance travel
- Tired, heavy, restless, aching legs



COMPRESSION LEVELS

- 15–20mmHg - Women’s 143, Men’s 183 (Three easy-to-shop sizes)
- 20–30mmHg - Women’s & Men’s 830 (Six measured sizes)



PINK
ARGYLE (45)



PURPLE
ARGYLE (46)



GRAPHITE
ARGYLE (47)



ROYAL BLUE
ARGYLE (48)



DARK NAVY
STRIPE (13)

FIBERS (*Latex-Free*)

- 87% Nylon
- 13% Spandex

MICROFIBER SHADES

SIZING

for women & men

830 | 20–30mmHg SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
SIZE	SS	SL	MS	ML	LS	LL

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 **WOMEN'S SHOE SIZES:** SHORT LEG 5.5–9.5, LONG LEG 8–12

Product Code Example: For product 832CMSW48	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color
	83	2	C	MS	W	48

143/183 | 15–20mmHg SIZING CHART

WOMEN'S SHOE	5–7	7.5–9.5	10–12
ANKLE CIRCUMFERENCE	7"–9" (18–23cm)	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)
MEN'S SHOE	6.5–8.5	9–11	11.5–14
ANKLE CIRCUMFERENCE	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)	8.5"–11.5" (21–29cm)
SIZE	A	B	C

Product Code Example: For women's product 143CB45	Series	Product Style	Size	Color
	143	C	B	45

Product Code Example: For men's product 183CB45	Series	Product Style	Size	Color
	183	C	B	45



SURGICAL AND MEDICAL SUPPLIES PTY LTD

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

Email: sales@surmed.com.au

Website: www.surmed.com.au

SIGVARIS is a registered trademark of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2017 SIGVARIS, Inc.

LIT-830SS 08.17