

# SURGICAL & MEDICAL SUPPLIES



Furniture Product Catalogue

**Surgical & Medical Supplies Pty Ltd** is a proudly South Australian owned and operated company supporting the Healthcare industry since 1970. Over the last 40 years Surgical & Medical Supplies has proven to be a leading supplier of medical equipment, distributing for some of the industries best known brands. Currently our customers can order from a wide variety of product ranges, including but not limited to:

Philips	Convatec/UnoMedical	Sigvaris	Select
Heine	Midmed	Dalcross	Alphacare


Initially established to service Theatre, Critical Care, and Hospital Supply Departments, Surgical & Medical Supplies has grown to provide a wide variety of Healthcare solutions. Our extensive product range now allows us to service customers across many disciplines. Some examples of this include:

Anaesthetics	Gynecology	Veterinary Practices
Podiatry	Urology	Day Surgery
General Practice	Emergency Medicine	Vascular
CSSD	Radiology	Scientific Laboratories
Respiratory Medicine	Dentistry	Aged care

Customers have the opportunity to save time and money by dealing directly with a supplier in Adelaide.

From our central locality on Kensington Road we offer a high level of service, training and support through:

- Extensive warehousing - customizable to suit your needs
- Local customer service
- Same day delivery options
- Immediate administrative support
- Sales and Marketing support
- Product and Clinical Training
- Demonstration and sample products
- Access to Product Specialists and Key Account Managers
- All just a short stroll from the bustling east end of Rundle Street.



Surgical & Medical Supplies is committed to providing International options and pricing, available with local service and support. All of this delivered to your door on the same day from the one locally owned business.

We appreciate your continued support.

**Surgical & Medical Supplies**  
**16 Kensington Road**  
**ROSE PARK SA 5067**

**8:30am—5pm**  
**Monday to Friday**



**Telephone: (08) 8332 1666**  
**Facsimile: (08) 8332 5132**  
**Website: [www.surmed.com.au](http://www.surmed.com.au)**  
**Email: [sales@surmed.com.au](mailto:sales@surmed.com.au)**

Proud supplier of leading brands



Convatec



**PHILIPS**



**HiCare™**



FS 604 754

# Contents

<b>PATIENT TROLLEYS</b>	Page	<b>RECLINERS / TRANSPORTERS</b>	Page
Contour Classic	6	Champion Recliners / Transporters	19
Contour Portare	6	Stephen H Anatomical/Flat Treatment	19
Contour Deluxe	7	Bariatric Heavy Care Patient Trolley	19
Contour Multi-X	7		
Contour Endosurge	8	<b>EMERGENCY &amp; UTILITY CARTS</b>	
Contour Ophthalmic	8	Aurion Emergency 3 Cart	20
Contour Endo-X	9	Aurion Operating Room Cart	20
		Aurion Utility Cart	20
<b>PATIENT TROLLEY ACCESSORIES</b>		Prolite Std Emergency Basic Body	21
Multi-X Accessories	10	Prolite Anaesthesia	21
Deluxe Accessories	10	Prolite Emergency	21
Standard Accessories	11-14	Prolite Procedure	21
<b>TRANSPORT &amp; TREATMENT CHAIRS</b>		<b>GENERAL WAITING ROOM FURNITURE</b>	22-23
Swifi Transfer Chair	15		
Deneo Blood Sampling Chair	15	<b>STAINLESS STEEL TROLLEYS</b>	
Ambeo Power Day Surgery Chair	15	Flat Top Trolley	24
		Two Shelf Trolley with Rails	24
<b>EXAMINATION COUCHES</b>		Single Drawer Trolley	24
Standard Examination Couch	16	Two Drawer Side by Side Trolley	24
Heavy Duty Examination Couch	16	Two Drawer Vertical Trolley	24
Mobile Examination Couch	16	Four Drawer Trolley	24
<b>GYNAECOLOGICAL COUCHES</b>		<b>EXAMINATION / CLINICAL ROOM EXTRAS</b>	
Three Piece Gynaecological Couch	17	Surgeon Stool	25
Two Piece Gynaecological Couch	17	Linen Skips	25
Elansa Gynaecological Couch	17	IV Stands	25
		X-Ray Viewer	26
<b>ELECTRIC EXAMINATION COUCHES</b>		Superlux LSS Halogen Light	26
All-Electric 2 Section	18	Magnifying Lamp	26
All-Electric 3 Section	18		
Physio 3 Section	18		

# Patient Trolleys

## Contour Classic

- Highly customisable for best patient care
- Multi-riser dropside options
- Height, width and length variations
- Lowest trolley height



Contour Classic Electric - Built for the patient

## Specifications

Length 2135mm  
 Length (Short) 1805mm  
 Width (Standard) 815mm  
 Width (Narrow) 705mm  
 Width (Wide) 965mm  
 Mattress (Standard) 1900 x 680mm  
 Mattress (Narrow) 1900 x 570mm  
 Mattress (Wide) 1900 x 830mm  
 Max. Load 300kg  
 Castor Size 200mm

### Electric

Standard Height 440 - 810mm

Floor Clearance for

Trolley Mover 545 - 915mm

Electric Backrest Optional

Electric Knee Break Optional

### Hydraulic

Extra Low Height 525 - 755mm

Low Height 560 - 825mm

High Height 625 - 955mm

Manual Knee Break Optional

## Contour Portare

- Highly manoeuvrable
- Aesthetic, easy to clean styling
- Entry-level transfer trolley
- Excellent low height



Contour Portare Electric - Effective patient transport

## Specifications

Length 2135mm  
 Width (Standard) 815mm  
 Width (Narrow) 705mm  
 Width (Wide) 965mm  
 Mattress (Standard) 1900 x 680mm  
 Mattress (Narrow) 1900 x 570mm  
 Mattress (Wide) 1900 x 830mm  
 Max. Load 300kg  
 Castor Size 200mm

### Electric

Standard Height 450 - 850mm

Electric Backrest Standard

Knee Break Optional

# Patient Trolleys

## Contour Deluxe

- Full length stainless steel dropsides
- Fold away for zero transfer gap
- Gas assisted for total ease of use
- Safety release at both ends



## Specifications

Length 2135mm  
 Width (Standard) 815mm  
 Width (Wide) 965mm  
 Mattress (Standard) 1900 x 680mm  
 Mattress (Wide) 1900 x 830mm  
 Max. Load 300kg  
 Castor Size 200mm

### Electric

Standard Height 580 - 955mm  
 Electric Backrest Optional  
 Electric Knee Break Optional

### Hydraulic

Standard Height 625 - 955mm  
 Manual Knee Break Optional

## Contour Multi-X

- Unobstructed head to foot tracking carriage
- Easily customised to suit all brands of X-Ray system
- Removable cassette holder for easy cleaning
- Fully translucent top for entire body coverage



## Specifications

Length 2135mm  
 Width (Standard) 815mm  
 Width (Narrow) 705mm  
 Width (Wide) 965mm  
 Mattress (Standard) 1900 x 680mm  
 Mattress (Narrow) 1900 x 570mm  
 Mattress (Wide) 1900 x 830mm  
 Max. Load 300kg  
 Castor Size 200mm

### Electric

Standard Height 545 - 915mm  
 Floor Clearance for  
 Trolley Mover 650 - 1020mm  
 Electric Backrest Optional  
 Electric Knee Break Not advised

### Hydraulic

Extra Low Height 615 - 845mm  
 Low Height 650 - 915mm  
 High Height Not advised  
 Manual Knee Break Optional

# Patient Trolleys

## Contour Endosurge

- Enhanced clearance base for camera
- Hydraulic or Electric operation
- Sliding top allows full body coverage
- Split dropside option



Contour Endosurge Electric - Increased patient coverage

## Specifications

Length 2135mm  
 Width (Standard) 815mm  
 Width (Narrow) 705mm  
 Width (Wide) 965mm  
 Mattress (Standard) 1900 x 680mm  
 Mattress (Narrow) 1900 x 570mm  
 Mattress (Wide) 1900 x 830mm  
 Max. Load 300kg  
 Castor Size 200mm  
 Sliding Top  
 Overall Movement 460mm  
 Centre to Head 175mm  
 Centre to Foot 285mm  
**Electric**  
 Standard Height 560 - 930mm  
 Electric Backrest NA  
 Electric Knee Break NA  
**Hydraulic**  
 Extra Low Height 525 - 755mm  
 Low Height 560 - 825mm  
 High Height 625 - 955mm  
 Manual Knee Break NA

## Contour Ophthalmic

- Infinitely adjustable headrest
- Touch button operation
- Excellent clearance for surgeon
- Scalloped head cushion



Contour Ophthalmic Electric - Easy patient positioning

## Specifications

Length 2210mm  
 Width (Standard) 815mm  
 Width (Narrow) 705mm  
 Width (Wide) 965mm  
 Mattress (Standard) 1900 x 680mm  
 Mattress (Narrow) 1900 x 570mm  
 Mattress (Wide) 1900 x 830mm  
 Max. Load 300kg  
 Castor Size 200mm  
**Electric**  
 Standard Height 455 - 825mm  
 Floor Clearance for  
 Trolley Mover 560 - 930mm  
 Electric Backrest Optional  
 Electric Knee Break Optional  
**Hydraulic**  
 Extra Low Height 525 - 755mm  
 Low Height 560 - 825mm  
 High Height 625 - 955mm  
 Manual Knee Break Optional



# Patient Trolleys

## Contour Endo-X

- Complete freedom of working height
- Extra distance between columns
- Split dropside and reversible chassis options
- Unimpeded camera movement



## Specifications

Length 2135mm  
Width (Standard) 815mm  
Width (Narrow) 705mm  
Width (Wide) 965mm  
Mattress (Standard) 1900 x 680mm  
Mattress (Narrow) 1900 x 570mm  
Mattress (Wide) 1900 x 830mm  
Max. Load 300kg  
Castor Size 200mm  
Electric  
Standard Height 560 - 930mm  
Electric Backrest NA  
Knee Break NA


# Patient Trolley Accessories

## Multi-X Accessories




**Laser Mark System** This is a line of numbers (6-0-6) bolted to the inside of both end rails of the trolley, coinciding with numbers (0-6) on both ends of the crossslide.

Laser lines from the X-Ray tube are then used to centre the cassette accurately.



**Ezi-Centa** - clips onto the cassette crossslide and by using the shadow cast from the tube light source it can be positioned quickly under target. The bar then rotates for unobstructed X-Ray.

**Buccì Support** A full length plastic support sheet can be fitted to allow other Buccì systems to fit and give total freedom of movement. A self-centering cassette holder is also available.



### Removable Chest X-Ray

A lightweight, all stainless frame Hooks onto the backrest to support any size of X-Ray cassette and grid. Height adjustable to assist with positioning.


### C-Spine Cassette Holder

This stainless frame holds all cassette sizes and drops into the IV sockets on the trolley. It is height adjustable with positioning along the top as well.

### Cutaway Corner

Two or four corners of the bucci can be modified to allow further travel of the cassette at the foot or head end.

## Deluxe Accessories



**X-Ray Subtop** A translucent subtop can be fitted to the Deluxe for full body X-Rays. Both ends of the subtop can be lifted for easy cleaning. A self-centring cassette holder is also available.

**Rail Clamp** Designed for the Deluxe, this clamp locks onto the side rail, near the backrest. It holds any accessory with a 5/8" spigot.

# Patient Trolley Accessories

## Standard Accessories



**Multi Riser Dropside & Reverse Folding** 6, 7 or 16 riser dropside can be fitted to some models with all dropside also able to be mounted folding from foot towards head.



**Flush Fold Dropside** Classic style dropside which folds down flush with trolley top, may be fitted, allowing slide transfer of mattress or horizontal X-Ray.



**Mid Position Catch** This catch locks the Classic style dropside at approx. 45 deg. To give more access to the patient and some continued support.



**Mattress** Three mattress thicknesses are available: 50mm, 75mm (standard issue) and 100mm. They are made from medium density foam and covered in pressure care material.



**100mm Pressure Care** This mattress is constructed of layered foam - slow recovery over medium density. It is covered with four way stretch pressure relieving material.



**Padded Bumpers** 25mm thick high density foam covered in the same material as the mattress. Fits both the Classic and Deluxe style dropside.



**Head End Accessory Sockets** A pair of sockets mount to the head end of a Contour top and hold accessories including end rails and Ophthalmic head board.



**Ophthalmic Head Board** Lightweight head board which drops into head end accessory sockets on a Contour Trolley. It tilts up and down and is lockable.



**Arm Board** The arm board is padded for comfort and is 450 x 200 with 250mm height adjustment. It also has a wrist break 120mm in from the end.



**Surgi Rail Short** The Surgi Rail Short is 400mm long.



**Surgi Rail Long** The Surgi Rail Long is 700mm long. It is constructed from 30 x 6 stainless and bolts to the underside of Contour tops.



**Surgi Rail Clamp** An alloy clamp which holds 5/8" bar and fits onto our Surgi rails.

# Patient Trolley Accessories

## Standard Accessories



**Knee Break** Wind handle or electric operated, the Break raises 150mm for patient comfort during recovery or long stay.



**Dropdown Push Handles** Strong, foldaway handles with rubber grip for pushing the trolley. 150mm high x 250mm between hand grips.



**Horizontal Oxy Holder** A fully stainless steel basket for 155 dia. x 500 long oxy bottles. Hung from head end of trolley and allows front loading of bottle.



**Drain Bag Holder** 270mm long stainless steel holder, fastened to the centre of a trolley top on both sides.



**Ventilator Support Brackets** A variety of brackets which mount to the head end of Contour tops to support oxy delivery systems.



**Patient Tie Point** Strong metal bracket mounted both sides of trolley top to attach customer tie strap.



**Easy Head Access** A cutaway head endrail can be fitted on to most Contour trolleys to allow greater access to the patients face and head.



**Hydraulic Height Increase** A 50mm spacer can be added to increase the trolley height and camera clearance.



**Document Holder** 380 x 40 x 225mm deep wire basket. Hangs from trolley end rails. May also be fitted to the dropdown foot board or folding monitor board.



**Footend Rail** 175mm high x 600mm wide stainless rail, drops into footend sockets for pushing from or used as patient feet stop.



**Foot Rail Tall** Fully welded stainless foot rail 235mm high and 600mm wide with one cross tube.



**Foot Rail Support** Fully welded stainless foot rail 345mm high and 600mm wide, with two cross tubes

# Patient Trolley Accessories

## Standard Accessories



**Foot Rail Oxy** Fully welded stainless steel foot rail with support rings and stop for a "C" size oxygen cylinder.



**Vertical Oxy Holder (Rail Mount)** Stainless steel basket for maximum bottle size of 120 dia. X mm deep.



**Vertical Oxy Holder (Corner Mount)** Stainless steel basket for maximum bottle size of 120 dia. X 320mm deep.



**Fixed Monitor Board** A 450 high x 600 wide all stainless fixed frame with a 20mm full surrounding lip for safety.



**Folding Monitor Board** A strong stainless frame with push handles. The board folds over and has two tie straps for fastening monitors 500 high x 650 wide.



**Dropdown Foot Board** Padded for patient comfort, the dropdown foot board has handles for pushing the trolley and folds down to increase the mattress length.



**Folding IV Pole** The folding I.V. can be mounted at both ends of most trolleys and generally folds down to the right. It can be installed to reverse fold.



**Round Pump Pole** 25 Dia. x 950mm fixed height stainless pole for securing pump units to. Bag holder welded to top.



**Meal Tray** A vinyl wrapped board scalloped for patient comfort, it can fit both Classic and Deluxe style dropsides 680 x 400 x 25 mm thick.



**Paper Roll Holder** A bracket is bolted onto the head end rail of the trolley with a lift out paper roll axel. The paper is for high patient turnover, minimizing sheet changing time.



**Ezi-Trak** Centre steering makes handling effortless whilst the wheels ride easily over bumps and hollow due to gas pressure.



**Tri-Control** Raise, lower or tilt your trolley from the foot end. Head end control can be achieved by reversing the base and using special decals.

# Patient Trolley Accessories

## Standard Accessories



**Anti-Static Strip** A special anti sparking rubber strip mounted under a base to discharge any buildup of static electricity.



**Removable Storage Tray** Plastic storage tray can be clipped on to the head or foot end of the base 385 x 300 x 130 mm deep.



**Base Oxy Holder** An all stainless fixed frame for a maximum bottle size of 140 dia x 890mm with a retaining strap.



**Electric Drive Compatibility** Enhanced clearance bases or custom push bars can be made to suit most electric drive units used today



**Tie Rail on Base** Solid stainless rail mounted on electric or end control only bases, for securing special patients.

# Transport & Treatment Chairs

## Swifi Transfer Chair



- 200kg capacity
- Adjustable arm rests
- Centralised brakes
- Retractable footrest
- Backrest & seat in class M1 flame retardant seamless covering
- Front shock absorbers

## Deneo Blood Sampling Chair



- 200kg capacity
- Fixed height & hydraulic variable height
- Light grey epoxy covered steel for the base & structure
- Two lateral integrated hand grips
- One piece structure seamless upholstery
- Two removable formed armrests with swivel adjusting system
- Seamless heat formed M1 covering

## Ambeo Power Day Surgery Chair



- Designed for every stage: Welcome, Preparation, Transfer, Procedure & Wake up
- Foam thickness 80mm
- 4 Electric motors: back & seatrest reclinable, variable height, foot cushion
- Hand control allows patient independence
- Built-in oxygen bottle holder
- Removable upholstery with antibacterial protection
- Weight capacity 200kg

# Examination Couches

## Standard Examination Couch



- Adjustable headrest
- Overall dimensions: 1820mmL x 610mmW x 780mmH
- Safe working load 125kg
- Vinyl upholstery, powder coated frame

## Heavy Duty Examination Couch



- Adjustable headrest
- Overall dimensions: 1820mmL x 610mmW x 780mmH
- Comfortable thick foam padding
- Safe working load 250kg
- Vinyl upholstery, powder coated frame

## Mobile Examination Couch



- Vinyl upholstery colours: grey or black
- Fixed 75mm thick foam padded top
- 125mm Rubber castors, 2 with locks
- Adjustable headrest
- Overall dimensions: 1820mmL x 610mmW x 780mmH
- Safe working load: 125kg
- Accessories available:
  - Safety rails
  - Wire basket
  - IV pole & clamp



## Gynaecological Couches

### Three Piece Gynaecological Couch



- Sturdy chrome frame
- Vinyl upholstery
- Gynae cut-out & folding footrest
- Adjustable headrest
- Overall dimensions: 1845mmL x 610mmW x 780mmH
- Safe working load 125kg

### Two Piece Gynaecological Couch



- Sturdy chrome frame
- Vinyl upholstery
- Gynae cut-out & folding footrest
- Adjustable headrest
- Overall dimensions: 1820mmL x 610mmW x 780mmH
- Safe working load 125kg

### Elansa Gynaecological Table



- Electric height adjustment from 46cm to 106cm
- Combined height & backrest movements
- Pedal control
- Sliding plastic pan
- Gas spring assisted backrest adjustments
- Bilateral receptacle for optional colposcope holder
- Seamless upholstery with class M1 flame retardant covering
- Removable seat
- Paper roll holder
- Weight Capacity 170kg

## Electric Examination Couches

### All-Electric 2 Section



- Overall dimensions: 1955mmL x 710mmW
- 200kg Safe working load
- Hand switch controlled electric motors for Hi-Lo lift function and backrest
- Optional dual control kit available - Footswitch for Hi-Lo lift, Hand switch for back
- Leveling feet for perfect stability
- Retractable wheels for easy mobility

### All-Electric 3 Section



- Overall dimensions: 1955mmL x 710mmW
- 200kg Safe working load
- 3 Motors for height, backrest and leg rest
- 6-Button hand switch controls all functions. Additional foot switch controls the height only
- Leveling feet for perfect stability
- Retractable wheels for easy mobility
- Super strong construction with heavier base frame to prevent tipping

### Physio 3 Section



- Unique surround bar activated electric motor for Hi-Lo function
- Gas strut assisted long section and short section
- Overall dimensions: 1955mmL x 710mmW
- 200kg Safe working load
- Strong, rigid frame with dual lifting points
- Quality upholstery finish in grey or blue vinyl
- Levelling feet for perfect stability
- Retractable wheels for easy mobility

## Recliners / Transporters

### Champion Recliners / Transporters



- Patients recline easily by themselves
- All steel, powder coated frame for durability
- Foldaway side table
- Push handle
- Retractable foot tray
- 5" Thermoplastic castors
- Quick release seat

### Stephen H Anatomical / Flat Hospital Treatment Chair



- Ideally suited for Chemotherapy & Renal Dialysis
- Headrest
- Painted stainless steel frame 304 grade with antimicrobial coat
- Motor for simultaneous movement of back rest & leg rest
- Movement by handset with easy access for the patient
- Trendelenburg 10°
- Cleanable fabric with antibacterial and antimycotic properties
- Flame retardant fabric
- Weight capacity 200kg for standard size
- Mid, Junior & Large also available

### Bariatric Heavy Care Patient Trolley



- Fully motor powered
- Designed specifically for overweight patients
- Electric assisted function with 4 programmed positions: Patient Trolley, Chair, Bed, CPR Function
- Motor driven 5th wheel for easy steering
- Lateral tilt of mattress base to improve patient mobility
- Safe working load 350kg

## Emergency & Utility trolleys

### Aurion Compact Emergency Carts

- Custom designed to suit all environments
- Lockable drawers with 3 central locking devices to choose from
- Individual tags for labelling each drawer
- Medication basin - stainless steel fully autoclavable
- Waste container
- 46 Accessories with personalized solutions
- Body composed of a single rounded piece of Baydur polyurethane
- Does not react with diluted acids, alkali and solvents
- Handles moulded into the body

### Aurion Emergency 3 Cart



### Aurion Operating Room Cart



### Aurion Utility Cart



# Emergency & Utility trolleys

## Prolite Trolleys

Designed around users to guarantee quickness of manoeuvres and time of intervention, safety, high hygiene standard, reliability, high level of customization and upgradeability during the years

### Prolite Standard Emergency Basic Body



- Overall measurements 835 x 645 x 940mmH
- Basic body includes: right side drugs / ampoules drawer
- Left side pullout shelf
- CPR cardiac massage board
- Rear oxygen bottle holder

### Prolite Anaesthesia



- Easy access to medications needed at the point of care (ABC's equipments)
- Security
- Organisation
- Ample workspace

### Prolite Emergency



- Emergency care to the patients as quickly as possible
- Easy manoeuvrability
- Simultaneous access
- Organisation

### Prolite Procedure



- Secure and portable storage option
- Keeps patient treatment supplies close at hand at all times
- Customizations
- Flexibility
- Ample work space

## General Waiting Room Furniture

We supply a range of seating to the Healthcare Industry

**Harmony Chrome**



**Soho**



**Harmony Square**



**Ventura Blue**



**Pixie Black**



**Patriot**



**Comtub**



**Safe 1a**



**Lindsay**



# General Waiting Room Furniture

**Vibe 2 Seater Chair with Arms**



**Pixie Linked**



**Harmony Beam**



**Harmony Beam 4**



**Vibe Chair Table Combo**



**Safe 2 Seater**



**Sofa Bed**



**Tillo Leather Grey**



# Stainless Steel Trolleys

- Made in Australia. All welding to ANZS-1554
- High quality 304 grade stainless steel
- Safety folds to all edges
- 100mm Rubber wheels
- Full extension telescopic ball bearing drawer runners

## Flat Top Trolley - no rails



50 x 50 x 90cmH  
60 x 50 x 90cmH



80 x 50 x 90cmH  
110 x 50 x 90cmH

## Two Shelf Trolley with rails



50 x 50 x 90cmH  
60 x 50 x 90cmH



80 x 50 x 90cmH  
110 x 50 x 90cmH

## Single Drawer Trolley



50 x 50 x 90cmH  
60 x 50 x 90cmH

## Two Drawer Side by Side Trolley



80 x 50 x 90cmH

## Two Drawer Vertical Trolley



50 x 50 x 90cmH  
60 x 50 x 90cmH

## Four Drawer Trolley



80 x 50 x 90cmH



## Examination / Clinical Room Extras

### Surgeon Stool



- Height adjustment 530mm to 720mm
- Hand control lever
- Strong and stable five star moulded base with 50mm castors
- 90mm thick high-density foam padded seat
- Seat diameter 380mm
- Vinyl over and under seat for easy cleaning
- Vinyl upholstery colours: black or grey
- Safe working load: 140kg
- Available with foot rest, back rest & foot control

### Linen Skips



- Powdercoated
- Tubular metal frame
- Cream powdercoat finish
- Overall dimensions: 460mmW x 540mmD x 870mmH
- Bag holder opening/internal dimensions: 395mmW x 355mmD
- 75mm Castors
- Double linen skip also available
- Stainless steel option

### IV Stands



- Full stainless steel construction or moulded plastic base
- Stable 5 star base with 2 hook pole
- Height: 1370mm to 2200mm
- 4 Prong IV poles available
- Base diameter: 500mmW
- Five 50mm twin wheel castors
- Stackable option

## Examination / Clinical Room Extras

### Single Bay X-Ray Viewer



- Overall dimensions: 377mmL x 105mmD x 500mmH
- Actual illuminated panel size: 335mmL x 400mmH
- 2 Key-hole slots in rear for wall mounting
- Can be sat upright on base, for use on desk or bench
- Supplied with plug lead that exits at bottom right of light box
- Excellent light diffusion and uniform illumination is achieved from: circular fluorescent tube, continuous opal Perspex viewing panel and internal white reflector surface
- Stainless steel strip contains grip rollers to hold X-rays
- Each individually switched bay is internally divided
- White powdercoated steel case is easy to clean
- Available in other sizes
- Available with LED lighting

### Superlux LSS Halogen Light




- 50 Watt 240 volt dichroic bulb
- Bulb included
- 1250 lux at one metre
- 44cm beam diameter
- Maximum reach 1070mm

### Magnifying Lamp



- 1.75X magnification through 120mm diameter diopter lens
- 22 watt circular fluorescent tube produces 1250 lumens of shadowless 'daylight' illumination
- Angle poise spring balanced adjustable arm. Full extension 1150mm
- Protective dust cover flap
- Comes complete with table clamp. Wall bracket or mobile stand optional





Surgical & Medical Supplies Pty Ltd  
16 Kensington Road  
Rose Park, SA 5067

Phone: 08 8332 1666

Fax: 08 8332 5132

Website: [www.surmed.com.au](http://www.surmed.com.au)

E-mail: [sales@surmed.com.au](mailto:sales@surmed.com.au)