


**We bring  
LIGHT  
into  
the dark.**



## Instrument Head for HEINE UniSpec® Tubes



Illumination by Fiber Optic cable and projector or 3.5V BETA Handle with illumination adaptor see accessories on page 113.

- ✦ **Annular Fiber Optic illumination.** Bright, reflex-free illumination.
- ✦ **Magnifying lens 1.5x.** High-resolution image.
- ✦ **Hinged viewing window with insufflation port.** Improved vision by insufflation.
- ✦ **Stainless Steel.**

### Instrument head for UniSpec Disposable Tubes

complete with viewing window and swivel lens 1,5x [01]	<b>E-003.18.100</b>
without viewing window and swivel lens	<b>E-003.18.098</b>
Viewing window only [02]	<b>E-003.18.099</b>
Swivel lens 1.5x, only (for sigmoidoscopes) [03]	<b>E-003.18.096</b>
same, but 1.8x (only suitable for proctoscopes and anoscopes)	<b>E-003.19.099</b>

## HEINE UniSpec® Disposable Tubes



85x20mm    130x20mm    250x20mm

- ✦ **Disposable.** No cross-infection.
- ✦ **No cleaning required.** Saves time.
- ✦ **Moulded in patient-friendly thermoplastic.**  
No cold sensation on insertion.
- ✦ **Light transmission through the tube wall.** Distal annular light exit.  
Bright, homogeneous illumination of the area examined, reflex reduction.
- ✦ **Graduated inside the tube for patient comfort.** Indication of insertion length.

### UniSpec Disposable Anoscopes (85 x 20 mm)

Pack of 25 disposable anoscopes	<b>E-003.19.925</b>
Box of 100 (4 packs of 25) anoscopes	<b>E-003.19.911</b>

### UniSpec Disposable Proctoscopes (130 x 20 mm)

Pack of 25 disposable proctoscopes	<b>E-003.19.825</b>
Box of 100 (4 packs of 25) proctoscopes	<b>E-003.19.811</b>

### UniSpec Disposable Sigmoidoscopes (250 x 20 mm)

Pack of 25 disposable sigmoidoscopes	<b>E-003.18.825</b>
Box of 100 (4 packs of 25) sigmoidoscopes	<b>E-003.18.811</b>

⊗: Single use only.

## HEINE Anoscope/Proctoscope Illumination Head



[01]

**One single illumination head for all the tubes shown on this page.**

For 2.5V or 3.5V handle or Fiber Optic illumination with Fiber Optic projector and cable.

- ✦ **A wide range of tubes with either direct or 90° direction of view available.** Many different diagnostic and treatment methods are possible.
- ✦ **Tubes can be rotated through 360°.** Select the ideal position for treatment or examination.
- ✦ **Viewing window and swivel lens both offer 2x magnification.** High-resolution image.
- ✦ **Smooth, rounded surface.** Easy to clean. F.O. version can be cleaned with conventional disinfectants.



[02]



[03]

Illumination head	2.5V XHL	3.5V XHL	F.O.
complete with swivel lens and viewing window [01]	<b>E-001.19.101</b>	<b>E-002.19.101</b>	<b>E-003.19.120</b>
Illumination head without accessories	<b>E-001.19.102</b>	<b>E-002.19.102</b>	<b>E-003.19.121</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.040</b>	<b>X-002.88.051</b>	
Swivel lens, only [02]			<b>E-000.19.103</b>
Viewing window, only [03]			<b>E-000.19.104</b>

## Tubes with straight direction of view

**Anoscope/Proctoscope with proximal illumination.**

- ✦ **Distal end of tube angled.** Assures a good view and easy access with instruments.
- ✦ **All-metal construction.** Easy to clean and disinfect with conventional solutions.

### Tubes with straight direction of view

Anoscope tube 60 mm x 22 mm dia., complete with obturator	<b>E-000.19.111</b>
Proctoscope tube 120 mm x 22 mm dia., complete with obturator	<b>E-000.19.107</b>



60x22 mm

120x22 mm

## Tubes with 90° direction of view (Blond)

**Proctoscope with proximal illumination.**

- ✦ **Polished metal viewing mirror.** Tough, unbreakable.
- ✦ **End-piece with mirror detachable.** Easy to clean, disinfect or sterilize.

### Proctoscope Tubes

Proctoscope tube, 90° direction of view, 130 mm x 25 mm dia.	<b>E-000.19.105</b>
Proctoscope tube, 90° direction of view, 100 mm x 20 mm dia.	<b>E-000.19.114</b>



100x20 mm

130x25 mm



Due to local regulations, some products may not be available in your country.



HEINE F.O. Anoscopes with annular Fiber Optic Illumination



70x25mm



70x20mm




70x15mm

Illumination by Fiber Optic cable and projector or 3.5V BETA handle with illumination adaptor see accessories on page 113.

- **Three different diameters: 15 mm, 20 mm and 25 mm.** Useful range of sizes.
- **Annular Fiber Optic illumination.** Bright illumination of the area to be examined.
- **Angled tube-end.** Ensures a good view and easy access with instruments.
- **Rounded proximal shield.** Easy insertion and rotation.
- **Tube can be dismantled.** Easy to clean and disinfect with conventional solutions.
- **Stainless steel.** Durable.

Anoscope	F.O.
70x25 mm dia., complete with obturator	E-003.19.325
70x20 mm dia., complete with obturator	E-003.19.320
70x15 mm dia., complete with obturator	E-003.19.315

## HEINE F.O. Proctoscopes with Distal Annular Illumination



Illumination by Fiber Optic cable and projector or 3.5V BETA handle with illumination adaptor see accessories on page 113.

- ✦ **Distal annular Fiber Optic (F.O.) illumination.** Bright, reflex-free illumination.
- ✦ **Magnifying lens 1.8x magnification.** High-resolution image.
- ✦ **Hinged viewing window with insufflation port.** Vision can be improved by insufflation.
- ✦ **Angled tube-end.** Improves vision and access.
- ✦ **End-piece can be unscrewed.** Easy to clean and disinfect with conventional solutions.
- ✦ **Stainless steel.** Durable.

Proctoscope	F.O.
<b>130x20mm dia.,</b> complete with obturator, viewing window, swivel lens and end-piece	<b>E-003.19.213</b>
Proctoscope tube without accessories [01]	<b>E-003.19.013</b>
Obturator, only [02]	<b>E-003.19.113</b>
End-piece, only [03]	<b>E-003.19.095</b>
Viewing window, only [04]	<b>E-003.18.099</b>
Swivel lens 1.8x, only [05]	<b>E-003.19.099</b>

## HEINE F.O. Sigmoidoscopes with Distal Annular Illumination

Illumination by Fiber Optic cable and projector or 3.5V BETA handle with illumination adaptor see accessories on page 113.



- ✦ **Available in various sizes: 300x20 mm, 250x20 mm, 200x20 mm, 250x16 mm, 200x11.4 mm.** Sizes suitable for most patients of all ages.
- ✦ **Distal annular Fiber Optic (F.O.) illumination.** Bright, reflex-free illumination.
- ✦ **Swivel lens with 1.5x magnification.** High-resolution image.
- ✦ **Opening viewing window with insufflation port.** Vision can be improved by insufflation.
- ✦ **Graduated in cm.** Indication of insertion length.
- ✦ **Plastic end-piece.** Insulated.
- ✦ **End-piece can be unscrewed.** Easy to clean and disinfect with conventional solution.
- ✦ **Stainless steel.** Durable.

Sigmoidoscope	F.O.
<b>300x20 mm dia.</b> with obturator, viewing window, swivel lens and end-piece [01]	<b>E-003.18.230</b>
End-piece, only	<b>E-003.18.095</b>
<b>250x20 mm dia.</b> with obturator, viewing window, swivel lens and end-piece	<b>E-003.18.225</b>
End-piece, only	<b>E-003.18.095</b>
<b>200x20 mm dia.</b> with obturator, viewing window, swivel lens and end-piece	<b>E-003.18.220</b>
End-piece, only	<b>E-003.18.095</b>
<b>250x16 mm dia.</b> with obturator, viewing window, swivel lens and end-piece [02]	<b>E-003.18.525</b>
End-piece, only	<b>E-003.18.088</b>
<b>200x11.4 mm dia.</b> with obturator, viewing window and swivel lens [03]	<b>E-003.18.240</b>
End-piece, only	<b>E-003.18.094</b>
Viewing window, only [04]	<b>E-003.18.099</b>
Swivel lens 1.5x, only [05]	<b>E-003.18.096</b>



## Double Bulb Insufflator/Hygiene Filter

For use with anoscope, proctoscope and rectoscope



- **Latex-free double bulb insufflator** for air insufflation.
- For possible use with anoscope, proctoscope and rectoscope.
- **Disposable hygiene filter** minimises contamination of the double bulb insufflator.

Double Bulb Insufflator	E-000.18.105
-------------------------	--------------

Disposable Hygiene Filters 25 pcs/pack	E-000.18.116
--	--------------

## Sponge Holder



- **High quality design in stainless steel.**

Sponge holder 400 mm working length	E-000.18.906
-------------------------------------	--------------

## HEINE Illumination Adaptor



Can be attached to any HEINE 3.5V power source to illuminate any HEINE fiber optic instrument, e.g. anoscope, proctoscope or sigmoidoscope.

<b>Illumination adaptor</b>	<b>3.5V XHL</b>
Illumination adaptor	X-002.98.002
XHL Xenon Halogen spare bulb	X-002.88.078



## Sigmoidoscope/Proctoscope Kit RE 7000



Set complete with: **Instrument head** for UniSpec Disposable Tubes with viewing window and swivel lens 1.5x  
 Pack of 25 UniSpec Sigmoidoscope Tubes 250 mm  
 Pack of 25 UniSpec Proctoscope Tubes 130 mm  
 Pack of 25 UniSpec Anoscope Tubes 85 mm  
 Sponge holder  
 Double bulb insufflator  
 F.O. Projector HK 7000, 150 W  
 Fiber Optic cable

<b>Sigmoidoscope/Proctoscope Kit RE 7000</b>	<b>F.O.</b>
Complete Sigmoidoscope/Proctoscope Kit RE 7000	<b>E-096.15.501</b>
Spare bulb, 150W	<b>Y-096.15.102</b>
Disposable Hygiene Filters 25 pcs/pack	<b>E-000.18.116</b>

230V is the standard mains voltage. If a different voltage is required, e. g. 120V, please specify when ordering.

## Sigmoidoscope/Proctoscope Kit BETA4 USB



Set complete with: **Instrument head** for UniSpec Disposable Tubes with viewing window and swivel lens 1.5x.  
 Pack of 25 UniSpec Sigmoidoscope Tubes 250 mm  
 Pack of 25 UniSpec Proctoscope Tubes 130 mm  
 Pack of 25 UniSpec Anoscope Tubes 85 mm  
 Sponge holder  
 Double bulb insufflator  
 Illumination adaptor  
 Connecting cord 2 m  
 BETA4 USB Rechargeable Handle with  
 E4-USB Medical approved plug-in power supply

<b>Sigmoidoscope/Proctoscope Kit BETA4 USB</b>	<b>3.5V XHL</b>
Complete Sigmoidoscope/Proctoscope Kit BETA4 USB	<b>E-097.27.388</b>
Spare bulb	<b>X-002.88.078</b>
Disposable Hygiene Filters 25 pcs/pack	<b>E-000.18.116</b>

