

Premium Stretcher Mattress

High performance mattress with outstanding durability and patient support for demanding acute environments.

Features

- 5-layered high performance core with up to a 300kg therapeutic support load
- Dual interface and central immersion regions with firm base support
- Integrated profile stabilizing sides and unique hinging system
- Double coated, welded polyurethane cover with fluid proof zip
- Superior 110mm mattress thickness
- Non-slip base material with high friction grip
- 5 year core foam and 4 year cover warranty
- Custom sizing and labeling available

Product Code





- Square Ends (electric stretcher); PT-CM-110-PRPF
- Truncated End (hydraulic stretcher); PT-CM-110-PRF

Mattress Size

- Standard; 680 W x 1900mm L
- Narrow; 570mm W x 1900mm L
- Wide; 830mm W x 1900mm L

Summary

The **Premia Forte** pressure relieving stretcher support surface is our high spec version of the standard Contour Premia mattress. The specialised core structure provides best protection to patients at risk of injury without compromising support in crucial areas. A Premiflex stretcher cover provides ultimate patient comfort and protection. This mattress delivers an effective outcome - combining both high risk pressure care and longevity for demanding stretcher applications.



17/04/2020