



Distributed in SA by



SURGICAL AND MEDICAL SUPPLIES PTY LTD

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

Email: sales@surmed.com.au

Website: www.surmed.com.au

Disposable Diathermy Electrodes



www.elitemedical.com.au

Diathermy Electrodes

The diathermy electrode range from Elite Medical consists of high quality single-use devices for a variety of surgical procedures.

They are designed and manufactured as a result of the demands of leading surgeons and specialists, ensuring clinical performance.

Non-stick PTFE coating (where supplied), reduces the build-up of carbonised material (eschar) on the electro-surgical tips during surgery. This reduces the need for higher output settings which can cause damage to the peripheral tissue.

CONTENTS




BLADE ELECTRODES

A range of non-stick PTFE blade electrodes available in standard and shielded styles.

Shielded tips prevent unwanted contact with tissue. They are designed to focus the energy on the target tissue, and offer a greater level of safety when working in cavities and areas of restricted vision.

Suit all common pencils, and supplied sterile, in boxes of 10.

Code	Description	Size
PMS2950	Blade Electrode Shielded Tip	63mm
PMS2920	Blade Electrode	63mm
PMS2951	Blade Electrode Shielded Tip	70mm
PMS2921	Blade Electrode	70mm
PMS2926	Blade Electrode Shielded Tip	100mm
PMS2925	Blade Electrode	100mm
PMS2928	Blade Electrode Shielded Tip	130mm
PMS2927	Blade Electrode	130mm
PMS2933	Blade Electrode Shielded Tip	165mm
PMS2929	Blade Electrode	165mm




NEEDLE ELECTRODES

A range of non-stick PTFE needle electrodes available in standard and shielded tip styles.

Shielded tips prevent unwanted contact with tissue. They are designed to focus the energy on the target tissue, and offer a greater level of safety when working in cavities and areas of restricted vision.

Suit all common pencils, and supplied sterile, in boxes of 10.

Code	Description	Size
PMS2972	Needle Electrode Shielded Tip	70mm
PMS2922	Needle Electrode	70mm
PMS2946	Needle Electrode Shielded Tip	100mm
PMS2945	Needle Electrode	100mm
PMS2935	Needle Electrode Shielded Tip	152mm
PMS2934	Needle Electrode	152mm



LOOP ELECTRODES

Designed with an ultra-fine tungsten wire for minimal eschar build up.

Suit all common pencils, and supplied sterile, in boxes of 5.

Code	Colour	Loop Size
PMS1742	Aqua	10mm x 10mm
PMS1746	Orange	15mm x 8mm
PMS1747	Pink	15mm x 12mm
PMS1743	Blue	15mm x 15mm
PMS1748	Green	20mm x 10mm
PMS1751	Grey	20mm x 12mm
PMS1745	Red	20mm x 15mm
PMS1744	Beige	20mm x 20mm




BALL ELECTRODES

Non-Stick PTFE coated ball electrodes minimise adhesion and prevent bleeders from resuming when concluding a procedure. This further allows for fulguration to be achieved with greater success rates than with standard stainless steel ball electrodes.

Suit all common pencils, and supplied sterile, in boxes of 10.

Code	Description	Ball Size	Length
PMS2905	Non-stick Ball	4mm	50mm
PMS2907	Non-stick Ball	4mm	110mm
PMS2907L	Non-stick Ball	4mm	130mm
PMS2906	Non-stick Ball	5mm	110mm
PMS2906L	Non-stick Ball	5mm	130mm



PENCIL ELECTRODES

All pencils come with pre-attached 3 metre cable.

Code	Description	UOM
PMS254	Single Use Hand Control Fingerswitch Pencil Button Switch Nonstick PTFE Coated Blade with Shielded Tip	Box/25
PMS252	Single Use Hand Control Fingerswitch Pencil Button Switch Nonstick PTFE Coated Blade	Box/25
PMS250	Single Use Hand Control Fingerswitch Pencil Button Switch Stainless Steel Blade	Box/25
240-014-3	Single Use Foot Control Pencil 8mm Switch Stainless Steel Blade	Box 10

Smoke Evacuation and Electrosurgical Generators

A range of associated operating theatre equipment is available for all electrosurgical applications, from large hospitals to small clinics.

For more information please contact us on 1800 354 836.





The perfect solution at your fingertips 24/7

Our extensive and trusted range of medical products will fulfil any requirement you may ever have while delivering improved patient outcomes, ease of use and superior value. We carefully source premium instruments from only the most respected manufacturers worldwide.

Electrosurgery

- Electrodes
- Pencils
- Bipolar Forceps
- Bipolar Cables & Adaptors
- Monopolar Forceps
- Monopolar Cables & Adaptors
- Insulated Speculums
- Generators & Smoke Evacuators
- Quivers & Holsters

Surgical Instruments

- Scalpels and Handles
- Scissors
- Forceps
- Suction Tubes
- Retractors & Skin Hooks
- Needle Holders
- Chisels and Gouges
- Fibre Optic Cables
- Laparoscopic Consumables

Instrument Care & Sterilisation

- Instrument Protection Caps
- Identification Tags
- Identification Tape
- Cleaning Brushes
- Cleaning Chemicals & Solutions

Patient Care & Safety

- Pressure Protection
- Operating Table Accessories
- Anaesthetics
- Workplace Health & Safety

About Us

Elite Medical is a family owned and operated company who, since October 1986, have been dedicated to serving the healthcare sector by providing high quality, cost effective and environmentally responsible surgical instruments, equipment and consumables.

Our Values

INTEGRITY - We act with honesty and transparency, always doing what is right by our stakeholders.

SERVICE - Our customers are our number one priority, therefore we respond to your enquiry with enthusiasm, professionalism and efficiency.

PROGRESS - We are continuously improving our processes in order to remain at the forefront of our industry.

Our People

We believe our people are the most important part of our business. Therefore our team is made up of the most passionate specialists in their field. Our business development managers are committed to searching out innovation and our customer service team go the extra mile to not only meet but exceed your expectations.

Our Community

Elite Medical is 100% Australian owned and operated so we have a practical understanding of your needs. In addition, we are dedicated to being good corporate citizens. In order to achieve this we regularly support charities such as Careflight and the Royal Flying Doctor Service. And we are always improving our practices to ensure they are environmentally responsible.

**Exactly the specialised medical supplies
you want, when you need them.**



SURGICAL AND MEDICAL SUPPLIES PTY LTD

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

Email: sales@surmed.com.au

Website: www.surmed.com.au