


Sterilisation Trays

Designed to protect delicate instruments during washing, sterilisation and storage. Constructed of a polymer resin to provide strength, durability and chemical resistance, whilst remaining lightweight and incredibly easy to use.

- Highly durable, heat-resistant plastic base and lid
- Silicone finger mat protects and helps organise instruments
- Slots in lid and base facilitate the circulation of steam and aid in the drying process
- Raised nodes keep silicone mat from sticking to the base
- Transparent lid for easy viewing of contents with snap-shut closure
- Sold as a full unit – base, lid, mat
- Compatible with steam and gas plasma sterilisation

Code	Dimensions
PST-1100SB	190.5mm x 63.5mm x 19mm
PST-2100SB	165mm x 102mm x 19mm
PST-3100SB	254mm x 152mm x 19mm
PST-3200SB	254mm x 152mm x 38mm
PST-3300SB	254mm x 152mm x 38mm (double layer with insert tray and mat)
FM-1810	462mm x 260mm




Silicone Finger Mat


Individual mat that can be cut to size for specific size needs. Versatile product may be used for replacements within other Sterilisation Trays.




FM-1810




Other variations available




PST-1100SB




PST-2100SB



PST-3100SB



PST-3200SB



PST-3300SB

Exactly the specialised medical supplies you want, when you need them.



SURGICAL AND MEDICAL SUPPLIES PTY LTD
 16 Kensington Road Rose Park SA 5067
 Telephone: (08) 8332 1666 Fax: (08) 8332 5132
 Email: sales@surmed.com.au
 Website: www.surmed.com.au