



Shower Bathing Trolleys



HAND-MADE WITH CARE
in Australia

MODSEL
nursing made easier

March 2021

Distributed in SA by



SURGICAL AND MEDICAL SUPPLIES PTY LTD

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

Email: sales@surmed.com.au

Website: www.surmed.com.au



Shower Bathing Trolleys

Our mobile shower trolleys are constructed for both patient safety and ease of use by allied health professionals, and feature a clean design with durable materials that are ideal for harsh environments.

Perfect for use as a shower bed or stretcher, these designs allow for a stable, reliable solution for easy bathing across the board. Some options in our range include designs specially made for bariatric patients, while others offer complete privacy features.



Aquarius Junior



Aquarius



Aquatuff Electric



Aquatuff Hydraulic



Aquatuff Bari



Custom Liner Colours



Change Table



Concealment Trolley

Aquarius Junior

A special junior sized shower trolley developed for small patients.

Features

- Economical hydraulic trolley with 180kg SWL
- Powder coated stainless steel dropsides and endrails
- 790 W x 1680mm Long overall size including bumpers
- Liner inner (min) base size - 560 W x 1450 L x 180mm Deep
- Adjustable end rails with 120mm extra length
- Standard hose length; 1100mm, custom lengths available
- Height Range - 680 to 880mm
- Quality 125mm twin wheel castors with total lock and steer
- 135mm lowest clearance under the base
- MADE IN AUSTRALIA

Product Code

- MST-SA-JUN

Summary

Working with children or small adults may not require a full sized shower trolley. This junior sized version is also easier to move around in tight facilities and for storage. The refined framework and proven mechanics awards long life value and consists of considered, easy to clean surfaces coupled with quality twin wheel castors to ensure the highest level of manoeuvrability.



Safe transfer



Easy drainage



Height adjustment



Length adjustment

Aquarius

A standard shower trolley designed to meet budget constraints.

Features

- Economical hydraulic trolley with 180kg SWL
- Powder coated stainless steel dropsides and endrails
- 795 W x 1950mm Long overall size including bumpers
- Liner inner (min) base size - 560 W x 1720 L x 180mm Deep
- Adjustable end rails with 120mm extra length
- Standard hose length; 1100mm, custom lengths available
- Height Range - 680 to 880mm
- Quality 125mm twin wheel castors with total lock and steer
- 135mm lowest clearance under the base
- MADE IN AUSTRALIA

Product Code

- MST-AQM

Summary

We are able to bring our Modsel shower trolley design quality to you at an affordable price, helping to meet tight budget constraints and maximising the number of working units available.

The refined framework and proven mechanics awards long life value and consists of considered, easy to clean surfaces coupled with quality twin wheel castors to ensure the highest level of manoeuvrability.



Safe transfer



Easy drainage



Height adjustment



Length adjustment

Aquatuff Electric

A stainless steel shower trolley with electrically operated 300kg SWL, wide top and backrest options.

Features

- Electric height adjustment with 300kg SWL
- 15 degree electric backrest option available
- Total stainless steel construction
- 795 W x 1950mm Long overall size including bumpers
- Liner inner (min) base size - 560 W x 1720 L x 180mm Deep
- Adjustable end rails with 120mm extra length
- 100mm wider comfort option (internal min width of 660mm)
- Standard hose length; 1100mm, custom lengths available
- Height Range - 680 to 900mm
- Quality 125mm twin wheel castors with total lock and steer
- MADE IN AUSTRALIA

Product Code

- Aquatuff Electric; MST-SA-EUO
- Electric backrest option - MST-AQTE-B
- 100mm wider top option - MST-SA-EWT

Summary

Having equipment that is robust and completely reliable is of the utmost importance when caring for a larger or heavier patients. We have over-engineered the componentry strength in this design and our decision to use total stainless steel reinforces the a lifetime of value that this bariatric shower trolley will give. The electric operation helps to safely lift 300kg a without overload, whilst four twin wheels with a brakes, electric backrest and wider top options all increase versatility in special needs.



Lies flat



Standard design



Safe transfer



Easy drainage



Height adjustment



Length adjustment

Aquatuff Hydraulic

A stainless steel shower trolley with hydraulically operated 300kg SWL and wide top option.

Features

- Hydraulic height adjustment with 300kg SWL
- Total stainless steel construction
- 795 W x 1950mm Long overall size including bumpers
- Liner inner (min) base size - 560 W x 1720 L x 180mm Deep
- Adjustable end rails with 120mm extra length
- 100mm wider comfort option (internal min width of 660mm)
- Standard hose length; 1100mm, custom lengths available
- Height Range - 660 to 900mm
- Quality 125mm twin wheel castors with total lock and steer
- MADE IN AUSTRALIA

Product Code

- Aquatuff Hydraulic; MST-AQT
- 100mm wider top option - MST-SA-EWT

Summary

Having equipment that is robust and completely reliable is of the utmost importance when caring for a larger or heavier patients. We have over-engineered the componentry strength in this design and our decision to use total stainless steel reinforces the a lifetime of value that this bariatric shower trolley will give. The hydraulic operation helps to safely lift a 300kg without overload, whilst four twin wheels with brakes and wider top all increase versatility in a special needs.



Safe transfer



Easy drainage



Height adjustment



Length adjustment



Wide top option

Aquatuff Bari

A robust, fully electric 300kg SWL shower trolley with extra wide base and top for Bariatric patients.

Features

- Electric height adjustment with 300kg SWL
- 15 degree electric backrest option available
- Total stainless steel construction
- 1195 W x 1950mm Long overall size including bumpers
- Liner inner (min) base size - 960 W x 1720 L x 180mm Deep
- Adjustable end rails with 120mm extra length
- Standard hose length; 1100mm, custom lengths available
- Height Range - 680 to 900mm
- Quality 125mm twin wheel castors with total lock and steer
- MADE IN AUSTRALIA

Product Code

- Aquatuff Bari; MST-SA-XWT
- Electric backrest option - MST-AQTE-B

Summary

Having equipment that is robust and completely reliable is of the utmost importance when caring for a larger or heavier patients. We have over-engineered the componentry strength in this design and our decision to use total stainless steel reinforces the a lifetime of value that this bariatric shower trolley will give. The electric operation helps to safely lift 300kg a without overload, whilst four twin wheels with a brakes with wider base and top all increasing a versatility in special needs.



Safe transfer



Easy drainage



Backrest option



Length adjustment

Custom Liner Colours

Customise the colour of your shower trolley liner from the standard Royal Blue.

Features

- Heavy duty washable, non-stick vinyl
- Ask for the available colours

Product Code

- MST-SA-CLC

Summary

There is a selection of light and dark colours available in the high grade vinyl used for the Aqua liners, that you can choose as an optional upgrade. Please contact us to discuss your requirements and view a colour chart.



Royal Blue



Vinyl colours

Change Table

A purpose built, non-slip change mat on a sturdy, height adjustable and safe patient platform.

Features

- 25mm high density foam core
- Heavy duty, welded vinyl cover
- 795 W x 1950mm Long overall size including bumpers
- Patient surface size of 600 W x 1820mm Long
- Both junior and wider sized options available
- Heights from 660 to 900mm (depending on trolley model)
- Hydraulic 180kg SWL to electric 300kg SWL options
- Quality 125mm twin wheel castors with total lock and steer
- Fits both Aquarius and Aquatuff trolley models
- MADE IN AUSTRALIA

Product Code

- MST-SA-CTU

Summary

This change mat expands the use of the Aquarius a trolley into a very safe, height adjustable change a table with multi-position, removable end rails and full length dropsides that fold back underneath the trolley top for ample work clearance. The mat a securely fastens underneath the trolley top to a prevent sliding and movement. This trolley is designed for harsh use and able to be rigorously cleaned.



Safe design



Customise

Concealment Trolley

A lightweight, removable frame and block-out, washable cover fitted to a trolley for transfer privacy.

Features


- Heavy duty washable vinyl cover
- Opens to full access either side and zips down each corner
- Smooth, one-piece, cleanable trolley top
- Patient surface size of 600 W x 1820mm Long
- 410mm clearance height from patient surface
- Heights from 660 to 900mm (depending on trolley model)
- Both junior and wider sized options available
- Hydraulic 180kg SWL to electric 300kg SWL options
- Quality 125mm twin wheel castors with total lock and steer
- Fits both Aquarius and Aquatuff trolley models
- MADE IN AUSTRALIA

Product Code

- MST-SA-CCU

Summary

Unzip one side and fold the cover back to slide a deceased patient onto the smooth, flat surface of a trolley top. The height adjustment and fully fold under dropsides help to ensure transfer is as easy as possible. The cover ensures total privacy and the mortuary trolleys are constructed for heavy washing and cleaning.



Fully enclosed



Customise

