

Stainless Steel Trolleys & Carts



CONTENTS

	Page
1. Instrument Trolleys	2
2. Dressing Trolleys	2
3. Powder Coated Trolleys	3
4. Procedure Cart	4
5. Resuscitation Carts	4
6. Anaesthetic Cart	4
7. Trolley Accessories	4

SA Distributor



SURGICAL AND MEDICAL SUPPLIES PTY LTD
16 Kensington Road Rose Park SA 5067
Telephone: (08) 8332 1666 Fax: (08) 8332 5132
Email: sales@surmed.com.au
Website: www.surmed.com.au

For all enquiries, contact
Chris Lee, Product Specialist

0476 777 767

clee@surmed.com.au

25/08/2021

Stainless Steel Trolleys - Instrument & Dressing

Common Features

- Fully welded
- Double safety folds to all shelf edges
- Ball bearing fork castors with grey rubber wheels
- 304 grade stainless steel construction - satin finish
- Full under shelf bracing support providing maximum stability and strength



Instrument Trolley Without Rails

Product Code	Dimensions W x D x H mm	Castor size mm
GD0870	490 x 490 x 900	100
GD0880	600 x 490 x 900	100
GD1040	750 x 490 x 900	100
GD0890	900 x 490 x 900	125
GD0920	1000 x 490 x 900	125
GD0930	1000 x 600 x 900	125
GD0940	1200 x 490 x 900	125
GD1050	1200 x 490 x 900	125



Instrument Trolley With Rails

Product Code	Dimensions W x D x H mm	Castor Size mm
GD0900	490 x 490 x 900	100
GD1060	600 x 490 x 900	100
GD0910	750 x 490 x 900	100
GD1070	900 x 490 x 900	125
GD1080	1000 x 490 x 900	125
GD1090	1000 x 600 x 970	125
GD0860	1200 x 490 x 900	125
GD1100	1200 x 600 x 970	125



Dressing Trolley with 1 Drawer
drawer height: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
GD0840	490 x 490 x 970	100
SST01LK	500 x 500 x 970 top locking drawer	125
GD0950	600 x 490 x 970	100
GD0850	750 x 490 x 970	100
GD0960	900 x 490 x 970	125



Dressing Trolley with 2 Drawers
drawer height: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
GD0970	490 x 490 x 970	100
SST02LK	500 x 500 x 970 top locking drawer	125
GD0370	600 x 490 x 970	100
GD0790	750 x 490 x 970	100
GD0980	900 x 490 x 970	125



Dressing Trolley with 3 Drawers
drawer height: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
GD0990	490 x 490 x 970	100
SST03LK	500 x 500 x 970 top locking drawer	125
GD1000	600 x 490 x 970	100
GD1010	750 x 490 x 970	100
GD0800	900 x 490 x 970	125



Dressing Trolley with 4 Drawers
drawer height: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
GD0810	490 x 490 x 970	100
GD0810CL	490 x 490 x 970 central locking	100
SST04LK	500 x 500 x 970 top locking drawer	125
GD0820	600 x 490 x 970	125
GD0830	750 x 490 x 970	125
GD1030	900 x 490 x 970	125



Dressing Trolley with 4 Drawers
drawer heights:
2x 125mm, 2x 250mm

Product Code	Dimensions W x D x H mm	Castor Size mm
SSAT04	600 x 490 x 970	125

Powder Coated Steel Trolleys - Instrument & Dressing

Common Features

- Fully welded
- Double safety folds to all shelf edges
- Ball bearing fork castors with grey rubber wheels
- White powder coated steel
- Full under shelf bracing support providing maximum stability and strength



Instrument Trolley With Rails

Product Code	Dimensions W x D x H mm	Castor Size mm
PCT00	500 x 500 x 970	95



Dressing Trolley with 1 Drawer
drawer height: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
PCT01	500 x 500 x 970	95



Dressing Trolley with 2 Drawers
drawer height: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
PCT02	500 x 500 x 970	95



Dressing Trolley with 3 Drawers
drawer height: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
PCT03	500 x 500 x 970	95



Dressing Trolley with 4 Drawers
drawer height: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
PCT04	500 x 500 x 970	95

Stainless Steel Carts

Common Features

- Fully welded
- Double safety folds to all shelf edges
- Ball bearing fork castors
- 304 grade stainless steel construction - satin finish
- Full under shelf bracing support providing maximum stability and strength



Procedure Cart With 6 Drawers

Drawer heights:
4x 150mm, 2x 200mm

Product Code	Dimensions W x D x H mm	Castor Size mm
250501	900 x 490 x 970	125



Anaesthetic Cart - 6 Drawers

1st & 4th Drawer - 75mm high
2nd & 3rd Drawer - 150mm high
5th & 6th Drawer - 100mm high

Product Code	Dimensions W x D x H mm	Castor Size mm
SSAT06	900 x 490 x 970	125



Resuscitation Cart with 2 Half Drawers

drawer heights: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
GC0420	800 x 500 x 940	100
GC3020	900 x 500 x 940	125



Resuscitation Cart with 4 Half Drawers

drawer heights: 125mm

Product Code	Dimensions W x D x H mm	Castor Size mm
GC0560	800 x 500 x 940	100
GC0600	900 x 500 x 940	125



Resuscitation Cart - 5 Drawers

4 half, 1 full width drawers, 125mm high

Product Code	Dimensions W x D x H mm	Castor Size mm
GC0580	750 x 490 x 970	125



Resuscitation Cart - 7 Drawers

6 half drawers, 125mm high
1 full width drawer, 250mm high

Product Code	Dimensions W x D x H mm	Castor Size mm
SSRT07	750 x 490 x 970	125

Trolley Accessories



**IV Pole with mounting clamp
GC2530**



**C-size Oxygen Cylinder Holder
EZ5660**

Can't find the size you want? We can have it custom made.

Contact Chris Lee, Product Specialist

0476 777 767 clee@surmed.com.au