

Distributed in SA by:



SURGICAL AND MEDICAL SUPPLIES PTY LTD

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

Email: sales@surmed.com.au

Website: www.surmed.com.au

Elite
MEDICAL
#because we care

Endoscopy Consumables



For enquiries, contact:
Chris Lee, Product Specialist
Surgical & Medical Supplies
m 0476 777 767
e clee@surmed.com.au

www.elitemedical.com.au

“We believe there is always a better way.”

For over 30 years, Elite Medical has stood for people and products of the highest quality to continually supply great value instruments and consumables for health service providers. We keep it simple giving the best quality, life-saving instruments, combined with responsive time-saving service.

While we are not the largest supplier or owned by a multi-national, we are an Australian family-owned business successfully providing the Australian Medical Industry. We know that we can make a difference, so we strive to achieve our partners' goals and dreams **#becausewecare**.

Why Elite Medical

With our focus on **continuous improvement**, we **look after our customers** seeking long-term partnerships, to conduct every transaction with **honesty and transparency**, **standing accountable** to our stakeholders and ourselves for best results.

Regardless of the size of your business or annual spend, our support will always be second to none, and you will always have a dedicated account manager to support you.

Simply order online, give us a call or join our subscription service. Speak to your account manager for more details.


Our Endoscopy Consumables range

Elite Medical's single-use endoscopy consumables offer you a broad range of high quality and cost-effective medical devices to meet and exceed your demands, reducing the risks of infections associated with reusable products.

Backed by over 30-years of experience in the market, our brands, **Kirurgia®** and **Neurein®**, provides innovative and reliable solutions engineered to adapt to existing practices, introducing a range of medical consumables designed to simplify clinicians' roles, tailored to the unique needs of Australians.



KIRURGIA
BY ELITE MEDICAL




NEUREIN®
INFECTION CONTROL SOLUTIONS BY ELITE MEDICAL



Customer Enquiries

1800 354 836

sales@elitemedical.com.au



Product Offerings

www.elitemedical.com.au

If you find an '' please contact us for special order products.*

Endoscopy Consumables

<i>Pre-Procedure</i>	<i>1</i>
- Disposable Endoscope Buttons and Valves	2
- Endoscope Water Bottle Tubing	4
- Irrigation Tubing	5
<i>Procedure</i>	<i>6</i>
- Retrieval Nets	7
- Snares	8
- Monopolar Cable	9
- Biopsy Forceps	10
- Injection Needles	11
- Bite Blocks	12
- Polyp Traps	13
- Xylocaine Pump Spray Nozzle	14
- Yankauer Handles Suction Tube	15
<i>Post-Procedure</i>	<i>16</i>
- Cleaning Brushes	17
- Protective Sponge	18
- Face Shields	19
- Surgical Gowns	20
<i>Special Order Products</i>	<i>21</i>



Pre-Procedure

Buttons, Valves and Tubing

Our products reduce the cost of CSSD reprocessing as well as eliminating the risk of cross-infection. All you need to get set up ready for your endoscopy list! We have a range of water bottle tubing, irrigation tubing and endoscope buttons and valves.

Disposable Endoscope Buttons and Valves

Buttons and valve kits together in a convenient package minimising waste. These are supplied sterile including tracking labels on the outer package.

- Simply open the package and insert them in your scope
- Supplied with tracking labels
- Single-use, sterile



Disposable Endoscope Buttons and Valves

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-BK-SAJ	Kirurgia® Endoscope Buttons and Valves Kit	Air, Water Suction & Biopsy Valve, With Water Jet, Compatible with Olympus® Endoscopes 140/145/160/165/180/185/190/240/260/290 Series, Single-use, Sterile	BOX/25
EMK-BK-SA*	Kirurgia® Endoscope Buttons and Valves Kit	Air, Water Suction & Biopsy Valve, Compatible with Olympus® Endoscopes 140/145/160/165/180/185/190/240/260/290 Series, Single-use, Sterile	BOX/25

Disposable Endoscope Buttons and Valves

Sterile single-use valves that solve the issue of reusable valve reprocessing, saving you time. Single-use valves ensure sterile valves for every procedure and eliminates the need to track valves with endoscopes.

- ⦿ Simply open the package and insert them in your scope
- ⦿ Supplied with tracking labels
- ⦿ Single-use, sterile



EMK-BKBV-B1



EMK-BK-NR

Disposable Endoscope Buttons and Valves

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-BKBV-B1	Kirurgia® Biopsy Valve	Compatible with Olympus® and Fujinon® Endoscopes, Single-Use, Sterile	BOX/100
EMK-BK-NR*	Kirurgia® Auxiliary Water Jet Connector	Compatible with Olympus® Endoscopes 140/160/ 180/190 Series, Single-Use, Sterile	BOX/100

Endoscope Water Bottle Tubing

Disposable water bottle and irrigation tubes supplied sterile to eliminate the risk of contamination and the cost of reprocessing.

- ⦿ Clamp on the water line to avoid leaking between procedures
- ⦿ Recommended to be used for up to 24 hours for standard procedures
- ⦿ Water bottle compatibility: Baxter®
- ⦿ Single use, sterile



EMK-BKS100165CO2U

Endoscope Water Bottle Tubing

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-BKS100145CO2U	Kirurgia® Endoscope Water Bottle Tubing	Compatible with Olympus® Endoscopes, Single-Use, Sterile	BOX/10
EMK-BKS100165CO2U	Kirurgia® Endoscope Water Bottle Tubing	Compatible with Fujinon® Endoscopes, Single-Use, Sterile	BOX/10
EMK-BKS100160U*	Kirurgia® Endoscope Water Bottle Tubing	Compatible with Pentax® Endoscopes, Single-Use, Sterile	BOX/10

Irrigation Tubing

Kirurgia® irrigation tubing is universally compatible with any GI Endoscope or sterile water bottle.



- ⦿ Rigid cap seals tightly for maximum efficiency
- ⦿ Eliminates the need for manual cleaning and reprocessing of reusable tubing
- ⦿ Water bottle compatibility: Baxter®
- ⦿ Recommended to be used for up to 24 hours for standard procedures
- ⦿ Single use, sterile

Irrigation Tubing			
CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-BK200230S	Kirurgia® Irrigation Tubing	3000mm LGTH, Universal Fitting, Compatible with Baxter® Water Bottles, Connects to any EGP-100 Irrigation Pump, Olympus® OFP, AFU-100®, and ERBE® EIP2 Pumps, Single Use, Sterile	BOX/10



Procedure

Instruments and Devices

We offer a selection of retrieval nets, polypectomy snares, biopsy forceps and injection needles suitable for gastroscopy and colonoscopy procedures.

These are offered in an ideal length [2300mm] for both upper and lower procedures. Offering a single size eliminates confusion in the procedure room and the risk of accidental wastage.

Retrieval Nets

Retrieval nets designed for easy retrieval of large polyps or foreign objects.

- ⦿ Standard length 2300mm
- ⦿ Complete range of net sizes 15-45 mm
- ⦿ Supplied with forceps
- ⦿ Supplied with tracking labels
- ⦿ Single-use, sterile



Retrieval Nets

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-40090006*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 25mm Opening Size, Oval Net, Single-use, Sterile	BOX/10
EMK-40090007*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 30mm Opening Size, Oval Net, Single-use, Sterile	BOX/10
EMK-40090065*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 35mm Opening Size, Octagonal Net, Single-use, Sterile	BOX/10
EMK-40090066*	Kirurgia® Rotatable Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 35mm Opening Size, Octagonal Net, Single-use, Sterile	BOX/10
EMK-40090097*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 15mm Opening Size, Oval Net, Single-use, Sterile	BOX/10
EMK-40090098*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 35mm Opening Size, Oval Net, Single-use, Sterile	BOX/10
EMK-40090099*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia X 2300mm LGTH, 45mm Opening Size, Oval Net, Single-use, Sterile	BOX/10

Snares

Hybrid hot snares suitable for gastroscopy and colonoscopy procedures.

- ⦿ Easy to use
- ⦿ Single-use, sterile



EMK-40020158

Snares			
CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-40020158	Kirurgia® Polypectomy Cold Snare	2.4mm dia. X 2300mm LGTH, 15mm Opening Size, Rhomb Shape, Single-use, Sterile	BOX/10
EMK-40020313*	Kirurgia® Hybrid [Hot] Polypectomy Snare	2.4mm dia. X 2300mm LGTH, 15mm Opening Size Oval Shape, Single-use, Sterile	BOX/10

Monopolar Cable

The Kirurgia® disposable monopolar cable suits most snares. This saves the expense of testing and reprocesing of reusable cables and ensures you always have a cable available.

- ⦿ 8mm Bovie Plug on machine end
- ⦿ 2.9mm socket on instrument end
- ⦿ 3m cable
- ⦿ Single use, sterile



EMK-220-084-3

Disposable Monopolar Cable

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-220-084-3	Kirurgia® Monopolar Cable	8mm Bovie Plug, Female, Single-use, Sterile	BOX/10

Biopsy Forceps

Single size range of biopsy forceps. Our range includes oval cup and alligator jaw, and are available with or without spikes.

- ⦿ Oval cup or alligator teeth
- ⦿ With & without spike
- ⦿ Suitable for both gastroscopy and colonoscopy procedures
- ⦿ Supplied with tracking labels
- ⦿ Single-use, sterile



ALLIGATOR JAW



EMK-40010070

OVAL CUP



EMK-40010027

Biopsy Forceps

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-40010027	Kirurgia® Endoscopy Biopsy Forceps	2.3mm dia. X 2300mm LGTH, Oval Cup, Non-Stick, w/o Spike, Single-use, Sterile	BOX/10
EMK-40010042*	Kirurgia® Endoscopy Biopsy Forceps	2.3mm dia. X 2300mm LGTH, Oval Cup, Non-Stick, w/ Spike, Single-use, Sterile	BOX/10
EMK-40010070*	Kirurgia® Endoscopy Biopsy Forceps	2.3mm dia. X 2300mm LGTH, Rat Tooth Alligator Jaw, w/o Spike, Single-use, Sterile	BOX/10
EMK-40010084*	Kirurgia® Endoscopy Biopsy Forceps	2.3mm dia. X 2300mm LGTH, Rat Tooth Alligator Jaw, w/ Spike, Single-use, Sterile	BOX/10

Injection Needles

Single size range of injection needles to eliminate confusion in the procedure room and the risk of accidental wastage.

- ⦿ 2300mm in length
- ⦿ Suitable for both gastroscopy and colonoscopy procedures
- ⦿ Supplied with tracking labels
- ⦿ Single-use, sterile



EMK-400140024
5mm

Injection Needles

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-40140024	Kirurgia® Endoscopy Injection Needle	2.4mm dia. X 2300mm LGTH, 25G X 5mm LGTH, Single-use, Sterile	BOX/10
EMK-40140025	Kirurgia® Endoscopy Injection Needle	2.4mm dia. X 2300mm LGTH, 25G X 6mm LGTH, Single-use, Sterile	BOX/10

Bite Blocks

Kirurgia® range of bite blocks come in adult and child sizes, either with or without straps to facilitate oxygenation.

- ⦿ Cost effective
- ⦿ Various sizes
- ⦿ Large pack quantities
- ⦿ Single-use



Bite Blocks

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-40040001	Kirurgia® Endoscopy Bite Block	60mm X 55mm, w/Oxygen Connector, Adult, Single-use	BOX/100
EMK-40040002	Kirurgia® Endoscopy Bite Block with Strap	60mm X 55mm, w/Oxygen Connector & Strap, Adult, Strap 500mm LGTH, Single-use	BOX/100
EMK-40040004	Kirurgia® Endoscopy Bite Block with Strap	85mm X 55mm, Adult, w/Blue Textile Strap, Strap 500mm LGTH, Single-use	BOX/100
EMK-40040012	Kirurgia® Endoscopy Bite Block	85mm X 45 mm, Child, Single-use	BOX/100
EMK-40040013	Kirurgia® Endoscopy Bite Block with Strap	85mm X 45 mm, w/Strap, Child, Strap 500mm LGTH, Single-use	BOX/100

Polyp Traps

Kirurgia® polyp traps have five clearly numbered collection chamber that lock into place.

- ⦿ Well designed lid
- ⦿ Easy to open
- ⦿ Five separate chambers in one removable receptacle
- ⦿ Supplied with tracking labels
- ⦿ Single-use, sterile



EMK-BKP-200

Polyp Traps

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-BKP-200	Kirurgia® Disposable Polyp Traps	55mm dia. X 230mm LGTH, Five Chambers, Single-use, Sterile	BOX/10

Xylocaine Pump Spray Nozzle

Kirurgia® single use Xylocaine pump spray nozzle comes in conveniently packed boxes of 50, they are supplied non sterile. These nozzles atomise the spray to provide good coverage.

- ⦿ Cost effective
- ⦿ Large pack quantities
- ⦿ Single-use



Xylocaine Pump Spray Nozzle

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-BK-QZ	Xylocaine Pump Spray Nozzle	Disposable Pump Spray Nozzle, Non-Sterile	BOX/100

Yankauer Handles Suction Tube

Individually packed Yankauer Handles supplied with crown tip, with or without vents.

- ⦿ Crown tip to prevent tissue damage
- ⦿ Available with or without vents for better suction control
- ⦿ Single-use, sterile



Yankauer Handles Suction Tube			
CODE	PRODUCT NAME	DESCRIPTION	QTY
AN170002	Mdevices® Yankauer Handle	8mm dia. Crown tip, w/o Vents, Single-use, Sterile	BOX/50
AN170003	Mdevices® Yankauer Handle	8mm dia. Crown tip, w/ Vents, Single-use, Sterile	BOX/50



Post-Procedure


Cleaning Room & Personal Protection Equipment

For your cleaning galley, we have a conveniently selected range of single-use endoscope cleaning brushes and sponges, along with personal protective equipment that includes single-use thumb looped gowns and long face visors.

Cleaning Brushes

Neurein® endoscope cleaning brushes are designed and supplied in various sizes. The combo brush kits are most popular for high-volume procedure lists.

- ⦿ Best value combo kit
- ⦿ Individually packed brushes
- ⦿ Can clean endoscopes post procedure
- ⦿ Single-use



Cleaning Brushes

CODE	PRODUCT NAME	DESCRIPTION	QTY
NR-EBF-1004	Neurein® Double-Ended Channel Cleaning Brush	6mm dia. Brush size X 2300mm LGTH, Single-use	BOX/50
NR-EBF-1011	Neurein® Double-Ended Combination Valve Channel Brush	11mm dia. / 6mm dia. Brush size X 2300mm LGTH, Single-use	BOX/100
NR-EBF-1024	Neurein® Double-Ended Channel Brush and Valve Brush Kit	6mm dia. / 6mm dia. Brush size X 2300mm LGTH, 11mm dia. / 5mm dia. Brush size X 160mm LGTH , Single-use	BOX/100
NR-EBF-1026	Neurein® Double-Ended Channel Brush and Valve Brush Kit	3mm dia. / 3mm dia. Brush size X 1200mm LGTH, 11mm dia. / 5mm dia. Brush size X 160mm LGTH, Single-use	BOX/50
NR-EBF-1031	Neurein® Double-Ended Valve Brush	11mm dia. / 5mm dia. Brush size X 160mm LGTH, Blue Plastic Handle, Single-use	BOX/200

Protective Sponge

Neurein® soft endoscope sponges are supplied individually packed and sterile, this enables the sponge to be dual use. It can be used to protect the distal tip prior to the procedure and also be used to clean the scope after the procedure.

- ⦿ Sponge has a dual use and can be used as a tip protector
- ⦿ Can clean endoscopes post procedure
- ⦿ Single-use, sterile



Protective Sponge			
CODE	PRODUCT NAME	DESCRIPTION	QTY
NR-BKA-SPA	Neurein® Sterile Endoscope Sponge Transportation & Cleaning	40mm dia. X 120mm LGTH, Single-use, Sterile	BOX/50

Face Shields

Neurein® medical face shields offer the wearer anti-fog direct splash protection.

- ⦿ Foam headband
- ⦿ Latex-free
- ⦿ Easily disposable



Face Shields			
CODE	PRODUCT NAME	DESCRIPTION	QTY
NR-BK1400XL	Neurein® Face Shields	XL 230mm LGTH, Single-use	BOX/100

Surgical Gowns

Neurein® washroom-disposable gowns have convenient thumb loops to stop the sleeves from rising, reducing the risk of skin contact with chemicals and bodily fluids.

- ⦿ Impermeable
- ⦿ Back ties
- ⦿ Thumb loops
- ⦿ Single Use



Surgical Gowns			
CODE	PRODUCT NAME	DESCRIPTION	QTY
NR-BKSG	Neurein® Barrier Gown w/Thumb Loops	XL, Single-use	BOX/100

Special Order Products

CODE	PRODUCT NAME	DESCRIPTION	QTY
EMK-BK-SA*	Kirurgia® Endoscope Buttons and Valves Kit	Air, Water Suction & Biopsy Valve, Compatible with Olympus® Endoscopes 140/145/160/165/180/185/190/240 /260/290 Series, Single-use, Sterile	BOX/25
EMK-BK-NR*	Kirurgia® Auxiliary Water Jet Connector	Compatible with Olympus® Endoscopes 140/160/180/190/ Series, Single-Use, Sterile	BOX/100
EMK-BKS100160U*	Kirurgia® Endoscope Water Bottle Tubing	Compatible with Pentax® Endoscopes, Single-Use, Sterile	BOX/10
EMK-40090006*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 25mm Opening Size, Oval Net, Single-use, Sterile	BOX/10
EMK-40090007*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 30mm Opening Size, Oval Net, Single-use, Sterile	BOX/10
EMK-40090065*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 35mm Opening Size, Octagonal Net, Single-use, Sterile	BOX/10
EMK-40090066*	Kirurgia® Rotatable Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 35mm Opening Size, Octagonal Net, Single-use, Sterile	BOX/10
EMK-40090097*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 15mm Opening Size, Oval Net, Single-use, Sterile	BOX/10
EMK-40090098*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia. X 2300mm LGTH, 35mm Opening Size, Oval Net, Single-use, Sterile	BOX/10
EMK-40090099*	Kirurgia® Endoscopy Retrieval Net	2.6mm dia X 2300mm LGTH, 45mm Opening Size, Oval Net, Single-use, Sterile	BOX/10
EMK-40020313*	Kirurgia® Hybrid [Hot] Polypectomy Snare	2.4mm dia. X 2300mm LGTH, 15mm Opening Size Oval Shape, Single-use, Sterile	BOX/10
EMK-40010042*	Kirurgia® Endoscopy Biopsy Forceps	2.3mm dia. X 2300mm LGTH, Oval Cup, Non-Stick, w/ Spike, Single-use, Sterile	BOX/10
EMK-40010070*	Kirurgia® Endoscopy Biopsy Forceps	2.3mm dia. X 2300mm LGTH, Rat Tooth Alligator Jaw, w/o Spike, Single-use, Sterile	BOX/10



The perfect solution for your business

Our extensive and trusted range of medical products will fulfil your requirements while delivering improved patient outcomes, ease of use and superior value. We carefully source from only the most respected manufacturers worldwide.

Surgical Consumables

- Diathermy Forceps & Vessel Sealers
- Diathermy Pencils, Tips & Accessories
- Smoke Evacuation Consumables
- Laparoscopic Consumables
- Endoscopy Consumables

Surgical Instruments

- Scalpels & Handles
- Scissors
- Forceps
- Suction Tubes
- Retractors & Skin Hooks
- Needle holders
- Chisels and Gouges
- Fibre Optic Cables

Instrument Care & Sterilisation

- Cleaning Brushes
- Instrument Protection Caps
- Tray Tags & Labels
- Baskets & Trays
- Wraps & Packaging
- Quality Monitoring & Testing
- Transportation & Disposal Aids

Patient Care & Surgical Safety

- Pressure Protection Aids
- Patient Positioning Pads
- Hand Immobilisers
- Steps Stools & Anti-Fatigue Devices
- Specialised Anaesthetics
- PPE

Elite
MEDICAL
#because we care

For any enquiry about the products on this catalogue or any of our extensive range, please contact us:

1800 354 836

sales@elitemedical.com.au

www.elitemedical.com.au

Contact the Day Surgeries Account Manager to discuss your requirements. We can visit your facility and work with you to provide the products you need. If there are products that your department require but you do not see them in this catalogue, please contact us as we have a wealth of knowledge for sourcing bespoke products and solutions.

Distributed in SA by:



SURGICAL AND MEDICAL SUPPLIES PTY LTD

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

Email: sales@surmed.com.au

Website: www.surmed.com.au

Elite
MEDICAL
#because we care

For any enquiry about the products on this catalogue or any of our extensive range, please contact us:

1800 354 836

sales@elitemedical.com.au

www.elitemedical.com.au