

SPECULUMS - REUSABLE INSULATED

Elite Medical offer an extensive range of insulated speculums ideally suited for LLETZ/LEEP procedures.

Sized for petite to obese patients, these specula increase patient comfort, enhance visualisation and security for the surgeon and offer a cost-effective solution.

All are insulated, autoclavable and have an integrated smoke evacuation tube.

Graves Insulated Speculums

Code	Description	Size
5206-07NC	Small	75mm x 20mm
5206-10NC	Medium	100mm x 40mm
5206-12NC	Large	120mm x 40mm
5206-15NC	Extra Large	150mm x 40mm
5206-18NC	Extra Extra Large	180mm x 44mm

Graves Zee-More Insulated Speculums

Zee-More Speculums are designed with a wider 40mm yoke area to provide a more exposed view of the cervix and easier access for use of additional instruments.

Code	Description	Size
5200-25NC	Small	90mm x 25mm
5200-30NC	Medium	100mm x 40mm
5200-35NC	Large	120mm x 40mm

Graves Zee-Max Insulated Speculums

Zee-Max Insulated Speculums are designed with a wider 70mm yoke area which is also rectangular to provide maximum visibility and instrument access.

Code	Description	Size
5200-40NC	Medium	100mm x 40mm
5200-45NC	Large	120mm x 40mm

Sims Double Ended Insulated Speculums

Code	Description	Size
F050-390	Small	55mm x 25mm / 75mm x 31mm
F050-395	Medium	175mm x 27mm / 85mm x 36mm
F050-400	Large	85mm x 39mm / 95mm x 43mm



Sims Double Ended Insulated Speculums. Also available in Single Use.


Distributed in SA by:




SURGICAL AND MEDICAL SUPPLIES PTY LTD
16 Kensington Road Rose Park SA 5067
Telephone: (08) 8332 1666 Fax: (08) 8332 5132
Email: sales@surmed.com.au
Website: www.surmed.com.au




For enquiries, contact:
Chris Lee, Product Specialist
m 0476 777 767
e clee@surmed.com.au



Graves Insulated Speculums



Graves Zee-More Insulated Speculums



Graves Zee-Max Insulated Speculums

