

Distributed in SA by:



**SURGICAL AND MEDICAL SUPPLIES PTY LTD**  
16 Kensington Road Rose Park SA 5067  
Telephone: (08) 8332 1666 Fax: (08) 8332 5132  
Email: [sales@surmed.com.au](mailto:sales@surmed.com.au)  
Website: [www.surmed.com.au](http://www.surmed.com.au)

**Elite**  
**MEDICAL**  
*#because we care*

*Disposable  
Diathermy Forceps*



For enquiries, contact:  
Chris Lee, Product Specialist  
Surgical & Medical Supplies  
m 0476 777 767  
e [clee@surmed.com.au](mailto:clee@surmed.com.au)

[www.elitemedical.com.au](http://www.elitemedical.com.au)



## ***“We believe there is always a better way.”***

For over 30 years, Elite Medical has stood for people and products of the highest quality to continually supply great value instruments and consumables for health service providers. We keep it simple giving the best quality, life-saving instruments, combined with responsive time-saving service.

While we are not the largest supplier or owned by a multi-national, we are an Australian family-owned business successfully providing the Australian Medical Industry. We know that we can make a difference, so we strive to achieve our partners' goals and dreams **#becausewecare**.

## ***Why Elite Medical***

With our focus on **continuous improvement**, we **look after our customers** seeking long-term partnerships, to conduct every transaction with **honesty and transparency**, **standing accountable** to our stakeholders and ourselves for best results.

Regardless of the size of your business or annual spend, our support will always be second to none, and you will always have a dedicated account manager to support you.

Simply order online, give us a call or join our subscription service. Speak to your account manager for more details.

## ***Our Sieger® Disposable Diathermy Forceps Range***


We've come to listen to the operation room professionals throughout the years, and we've found manufacturers that will deliver the highest quality of disposable forces for your specific needs. Regardless of your surgery specialty, Elite Medical's extensive and trusted range of Disposable Diathermy Forceps will fulfil your electrosurgical needs.



*Customer Enquiries*

1800 354 836

sales@elitemedical.com.au



*Product Offerings*

[www.elitemedical.com.au](http://www.elitemedical.com.au)

## ***Disposable Diathermy Forceps***

<i>Our Brands</i>	5
<i>Sieger Blue®</i>	6
<i>Sieger Blue® Bipolar Forceps</i>	7
<i>Sieger Blue® Bipolar Forceps - No Cable</i>	9
<i>Sieger Blue® Monopolar Forceps</i>	10
<i>Sieger Blue® Adaptors and Cables</i>	11
<i>Sieger Non-Stick® Forceps</i>	12
<i>Gen2® Forceps</i>	14
<i>Sieger Platinum® Forceps</i>	16



# SIEGER

BY ELITE MEDICAL

With our over 30 years of experience providing service to the Australian health services, we've meticulously curated a selection of Australia's most demanded forceps patterns.

We've harnessed trust amongst medical staff across the country, and we have heard their demands. Because we care, Sieger® Diathermy Forceps product lines result from years of development to provide the best solutions to electrosurgery professionals to provide the greatest value with the highest quality of materials for your needs.

# Our Brands



**SIEGER**  
BLUE



**SIEGER**  
PLATINUM



**SIEGER**  
NON-STICK



**Gen<sup>2</sup>**  
BY SIEGER

# SIEGER BLUE

## *The cost effective, safe alternative*

The Sieger Blue® range of forceps is the original design and is available in many patterns and sizes, including bipolar and monopolar options. With over 34 configurations to cover most requirements, you can depend on Sieger Blue® to meet your needs.

Our extensive range of monopolar and bipolar forceps are crafted with the highest quality, providing a seamless experience as reusable instruments, conveniently removing the need for cleaning and reducing the chance of contamination. Available with or without pre-attached cable.

### Tips Sizes

0.5  
mm



1.0  
mm



1.5  
mm



2.0  
mm



Smooth Tip



Non-Stick Tip



Bipolar Plug



Serrated Tip



Toothed [1x2] Tip



Monopolar Plug

# Bipolar Forceps



Sieger Blue®  
McPherson Bipolar Forceps

CODE	DESCRIPTION	QTY
100-005-09050	0.5mm Smooth Straight Tip X 90mm LGTH, Single-use, Sterile	BOX/10



Sieger Blue®  
Iris Bipolar Forceps

CODE	DESCRIPTION	QTY
100-045-11050	0.5mm Smooth Straight Tip X 100mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



Sieger Blue®  
Jewellers Bipolar Forceps

CODE	DESCRIPTION	QTY
100-085-11030	0.3mm Smooth Tip X 114mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-085-11050	0.5mm Smooth Tip X 114mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



Sieger Blue®  
Adson Bipolar Forceps

CODE	DESCRIPTION	QTY
100-115-13050	0.5mm Smooth Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-115-13050S	0.5mm Serrated Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-115-13100	1.0 mm Smooth Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-115-13100S	1.0 mm Serrated Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/10

# Bipolar Forceps



**Sieger Blue®  
Cushing Bipolar Forceps**

CODE	DESCRIPTION	QTY
100-135-18100	1.0mm Smooth Tip X 180mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-135-18100S	1.0mm Serrated Tip X 180mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



**Sieger Blue®  
Jansen Bayonet Bipolar Forceps**

CODE	DESCRIPTION	QTY
100-185-20100	1.0mm Smooth Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-185-20100S	1.0 mm Serrated Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-185-20200	2.0 mm Smooth Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-185-20200S	2.0mm Serrated Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-185-22100	1.0mm Smooth Tip X 220mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-185-22100S	1.0mm Serrated Tip X 220mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-185-22200S	2.0mm Serrated Tip X 220mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



**Sieger Blue®  
Jansen Bayonet Bipolar Forceps  
- Angled Shaft**

CODE	DESCRIPTION	QTY
100-195-20100	1.0mm Smooth Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-195-20100S	1.0mm Serrated Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



**Sieger Blue®  
Jansen Bayonet Bipolar Forceps  
- Angled Tip**

CODE	DESCRIPTION	QTY
100-196-22150	1.5mm Smooth Tip X 216mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



# Bipolar Forceps - No Cable



Sieger Blue®  
Jewellers Bipolar Forceps  
No Cable

CODE	DESCRIPTION	QTY
100-084-11030	0.3mm Smooth Tip X 114mm LGTH, Single-use, Sterile	BOX/10
100-084-11050	0.5mm Smooth Tip X 114mm LGTH, Single-use, Sterile	BOX/10



Sieger Blue®  
Adson Bipolar Forceps  
No Cable

CODE	DESCRIPTION	QTY
100-114-13100	1.0mm Smooth Tip X 127mm LGTH, Single-use, Sterile	BOX/10
100-114-13100S	1.0mm Serrated Tip X 127mm LGTH, Single-use, Sterile	BOX/10



Sieger Blue®  
Straight Bipolar Forceps  
No Cable

CODE	DESCRIPTION	QTY
100-134-17100	1.0mm Smooth Straight Tip X 170mm LGTH, Single-use, Sterile	BOX/10



Sieger Blue®  
Bayonet Bipolar Forceps  
No Cable

CODE	DESCRIPTION	QTY
100-184-20100	1.0mm Smooth Tip X 200mm LGTH, Single-use, Sterile	BOX/10

# Monopolar Forceps



Sieger Blue®  
Adson Monopolar Forceps

CODE	DESCRIPTION	QTY
150-036-13100	1.0mm Smooth Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
150-036-13100S	1.0mm Serrated Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



Sieger Blue®  
Potts-Smith Monopolar Forceps

CODE	DESCRIPTION	QTY
150-046-18100S	1.0mm Serrated Tip X 180mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
150-046-20100S	1.0mm Serrated Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
150-046-24100K	1.0mm Toothed Tip X 240mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
150-046-24100S	1.0mm Serrated Tip X 240mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
150-046-24200S	2.0mm Serrated Tip X 240mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



Sieger Blue®  
Jansen Bayonet Monopolar Forceps

CODE	DESCRIPTION	QTY
150-056-19150S	1.5mm Serrated Tip X 190mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
150-056-22150S	1.5mm Serrated Tip X 220mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



Sieger Blue®  
Jansen Bayonet Monopolar Forceps  
Angled Tip

CODE	DESCRIPTION	QTY
150-066-22200S	2.0mm Serrated Tip X 220mm LGTH, w/ Cable, Single-use, Sterile	BOX/10

# Adaptors and Cables



Sieger Blue®  
Monopolar Adaptor

CODE	DESCRIPTION	QTY
ALC-634	4mm banana plug to 8mm Bovie	EACH



Sieger Blue®  
Monopolar Single Use

CODE	DESCRIPTION	QTY
220-024-3ES	8mm Bovie Pin, 3m LGTH	BOX/10



Sieger Blue®  
Bipolar Cable



Sieger Blue®  
Bipolar Cable

CODE	DESCRIPTION	QTY
190-014-3	3m LGTH, US 2 Pin Round Angled Shrouded plug, Single-use	BOX/10
190-014-3ST	3m LGTH, US 2 Pin Angled Non-Shrouded plug, Single-use	BOX/10

# SIEGER

## NON-STICK

### *The choice for Neuro, E.N.T. and Spinal Surgery*

Our unique bio-compatible silver alloy tip coating diverts heat quickly from the tips, preventing localised tissue adhesion, requiring less cleaning during surgery and reducing the overall amount of toxic surgical smoke. With its distinctive black and yellow colouring, the **Sieger Non-Stick®** range consists of a range of forceps with silver alloy coating, designed to minimize eschar build-up and provide a non-stick surface. Supplied with pre-attached cable.

#### Standard Forceps



'Hot spots' occur during coagulation



'Hot spots' cause localised sticking



Tissue seared to the tip 'breaks off'

#### Non-Stick Forceps



'Hot spots' start to occur during coagulation



Heat is rapidly conducted away to prevent adhesion





Forceps removed without adhesion

#### Non-Stick Tip



#### Tips Sizes



#### Fixed Pin Plug



# Non-Stick Forceps



**Sieger Non-Stick®  
Adson Bipolar Forceps**

CODE	DESCRIPTION	QTY
100-115-13050NS	0.5mm Non-Stick Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-115-13100NS	1.0mm Non-Stick Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



**Sieger Non-Stick®  
Bayonet Bipolar Forceps**

CODE	DESCRIPTION	QTY
100-185-20050NS	0.5mm Non-Stick Straight Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-185-20100NS	1.0mm Non-Stick Straight Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
100-185-22100NS	1.0mm Non-Stick Straight Tip X 220mm LGTH, w/ Cable, Single-use, Sterile	BOX/10



**Sieger Non-Stick®  
Angled-Shaft Bayonet Bipolar Forceps**

CODE	DESCRIPTION	QTY
100-195-17050NS	0.5mm Non-Stick Straight Tip X 170mm LGTH, w/ Cable, Single-use, Sterile	BOX/10

# Gen<sup>2</sup>

BY SIEGER

## Eco-friendly Forceps


The innovative **Gen2**<sup>®</sup> product range combines smart design with advanced manufacturing techniques. With a focus on sustainability and versatility, this range of forceps represents the latest in diathermy forceps design while offering sustainability and outstanding value. Supplied with pre-attached cable.



Smooth Tip



Serrated Tip



Tip Sizes



Flexi Pin Plug



# Sustainable Innovative Forceps



**Gen2®  
Adson Bipolar Forceps**

CODE	DESCRIPTION	QTY
3121-05	0.5mm Smooth Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/25
3121-05S	0.5mm Serrated Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/25
3121-10	1.0mm Smooth Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/25
3121-10S	1.0mm Serrated Tip X 127mm LGTH, w/ Cable, Single-use, Sterile	BOX/25



**Gen2®  
Jewellers Bipolar Forceps**

CODE	DESCRIPTION	QTY
2111-05	0.5mm Smooth Tip X 114mm LGTH, w/ Cable, Single-use, Sterile	BOX/25



**Gen2® Angled-Tip  
Jewellers Bipolar Forceps**

CODE	DESCRIPTION	QTY
2111-05A	0.5mm Smooth Angled Tip X 114mm LGTH, w/ Cable, Single-use, Sterile	BOX/25

# SIEGER PLATINUM

## Fully disposable forceps with superior non-stick properties

The Sieger Platinum range of forceps have been specifically designed for delicate neurosurgical procedures that require extra precision, performance, and balance.

Manufactured to the highest specifications, you can depend on Sieger Platinum!

Sieger Platinum® disposable forceps feature highly effective gold coated silver alloy tips for outstanding non-stick properties and are of such high-quality they feel the same as reusable instruments. Supplied with pre-attached cable. Ready to use with no need to clean, sterilise or check insulation.

### Standard Forceps



'Hot spots' occur during coagulation



'Hot spots' cause localised sticking



Tissue seared to the tip 'breaks off'

### Gold Coated Silver Alloy Tip Forceps



'Hot spots' start to occur during coagulation



Heat is rapidly conducted away to prevent adhesion




Forceps removed without adhesion

### Gold Coated Silver Alloy Tip



Tip Sizes



Flexi Pin Plug





# Gold Coated Silver Alloy Forceps



## Sieger Platinum® Bayonet Bipolar Forceps

CODE	DESCRIPTION	QTY
9201-10NS	1.0mm Non-Stick with Silver Alloy Tip X 200mm LGTH, w/ Cable, Single-use, Sterile	BOX/10
9221-10NS	1.0mm Non-Stick with Silver Alloy Tip X 220mm LGTH, w/ Cable, Single-use, Sterile	BOX/10

Distributed in SA by:



**SURGICAL AND MEDICAL SUPPLIES PTY LTD**

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

Email: [sales@surmed.com.au](mailto:sales@surmed.com.au)

Website: [www.surmed.com.au](http://www.surmed.com.au)

**Elite**  
**MEDICAL**  
*#because we care*

*For any enquiry about the products on this catalogue or  
any of our extensive range, please contact us:*

**1800 354 836**

**[sales@elitemedical.com.au](mailto:sales@elitemedical.com.au)**

**[www.elitemedical.com.au](http://www.elitemedical.com.au)**