


**AQ Medical Vaccine Fridges** are a **robust, cost effective** solution for all of your **vaccine refrigeration** needs.



### Features

- Self closing door to ensure temp never breaches vaccine temps
- Temperature keeps from 2°C - 8°C accurately and evenly
- High-low temp audible and visual alarm
- Suitable for small vaccine bottles
- Factory installed lock ensures security
- Electronic temperature controller and display
- 24 months warranty with national service support



## AQM 280

<b>Dimensions (HxWxD)</b>	<b>1625 x 580 x 584 mm</b>
<b>Capacity</b>	<b>280L</b>
<b>Temp Range</b>	<b>2°C - 8°C</b>
<b>No. of Shelves</b>	<b>6</b>
<b>Warranty (parts and labour)</b>	<b>24 Months</b>
<b>Power Supply</b>	<b>10amp</b>
<b>Feet Included</b>	<b>Feet at front, rollers at back</b>

SA Distributor



**SURGICAL AND MEDICAL SUPPLIES PTY LTD**  
16 Kensington Road Rose Park SA 5067  
Telephone: (08) 8332 1666 Fax: (08) 8332 5132  
Email: [sales@surmed.com.au](mailto:sales@surmed.com.au)  
Website: [www.surmed.com.au](http://www.surmed.com.au)

 **1800 658 120**

 **[avemquirks.com.au](http://avemquirks.com.au)**