


AQ Medical Vaccine Fridges are a **robust, cost effective** solution for all of your **vaccine refrigeration** needs.



Features

- Self closing door to ensure temp never breaches vaccine temps
- Temperature keeps from 2°C - 8°C accurately and evenly
- High-low temp audible and visual alarm
- Suitable for small vaccine bottles
- Factory installed lock ensures security
- Electronic temperature controller and display
- 24 months warranty with national service support



AQM 350

Dimensions (HxWxD)	1900 x 580 x 584 mm
Capacity	350L
Temp Range	2°C - 8°C
No. of Shelves	6
Warranty (parts and labour)	24 Months
Power Supply	10amp
Feet Included	Feet at front, rollers at back

SA Distributor



SURGICAL AND MEDICAL SUPPLIES PTY LTD
16 Kensington Road Rose Park SA 5067
Telephone: (08) 8332 1666 Fax: (08) 8332 5132
Email: sales@surmed.com.au
Website: www.surmed.com.au

 **1800 658 120**

 **avemquirks.com.au**