


Part No. ML201-LED

Standard LED Maggylamp Data Sheet



Fitted with	LED PCB
CRI	>80
Light output	1000 lumen
Colour temperature	5000K
Lens	Biconvex glass, Ø 127 mm
Magnification	2 x, 4 dioptre
Lamp shade	ABS plastic
Arms	chrome flex arms, 380mm long
Base housing	powder coated cast iron
Power lead	2m, AUS/NZ plug
Connected load	240V AC 50/60 Hz
Power consumption	approx. 12 Watt
Weight net	6.5kg
Weight carton	7kg
Package dimensions	L47 x W27 x H34cm

DIMENSIONS (Millimetres)



Data may vary due to component tolerances.
E & OE

PH:1300 788 239
www.redbank.net.au
©07.2017 Redbank Group


SURGICAL AND MEDICAL SUPPLIES PTY LTD
16 Kensington Road Rose Park SA 5067
Telephone: (08) 8332 1666 Fax: (08) 8332 5132
Email: sales@surmed.com.au
Website: www.surmed.com.au


lighting