



# Medisafe Plus 1000



**Medisafe Plus 1000** is one of the most feature-packed models of vaccine fridges available in Australia.

Medisafe Plus 1000 medical fridge is a purpose-built vaccine fridge designed and built in Australia in line with the National Vaccine Storage guidelines "Strive for 5".

Medisafe Plus vaccine fridges are feature packed, resulting in a vaccine storage solution that gives you peace of mind – knowing your vaccines are stored in the safest possible manner.

The Medisafe Plus 1000 is QCPP accredited and TGA approved.

The cloud monitoring acts as an automated data logger which allows you to view and print your temperature data in real time, from any location, on any web enabled smartphone, tablet, or computer.

All units also the ability to have a BMS system plugged into the back of the unit for hospital remote monitoring.

## Features

- 24 months Parts and Labour warranty
- Dedicated 24-hour technical support
- QCPP cold chain accredited
- TGA approved
- Min/Max Controller
- Temperature Controller
- Dual temperature probes for ambient and simulated vaccine temperatures
- AQ Box Cloud Monitoring – with SMS and email alerts of power or temperature failure
- Audio/Visual alarms
- Lock and Key
- Adjustable shelves
- Low Energy Consumption

## Optional Extras

- UPS Unit with Back up battery
- Preventative Maintenance program with extended warranty

<b>DIMENSIONS (HxWxD)</b>	<b>2200 x 1130 x 745</b>
<b>GROSS WEIGHT</b>	<b>210 Kg</b>
<b>GROSS CAPACITY</b>	<b>980 Lt.</b>
<b>TEMPERATURE RANGE</b>	<b>2° to 8°C</b>
<b>NUMBER OF SHELVES</b>	<b>10</b>
<b>POWER CONSUMPTION PER 24Hr</b>	<b>Tec = 3.57 kWh MEPS = 3.29 kWh</b>
<b>REFRIGERANT</b>	<b>R290a</b>

Local distributor:



**SURGICAL AND MEDICAL SUPPLIES PTY LTD**  
16 Kensington Road Rose Park SA 5067  
Telephone: (08) 8332 1666 Fax: (08) 8332 5132  
Email: [sales@surmed.com.au](mailto:sales@surmed.com.au)  
Website: [www.surmed.com.au](http://www.surmed.com.au)