



Medisafe Plus 80



Medisafe Plus 80 is the smallest and most convenient unit in the Medisafe Plus range.

Medisafe Plus 80 medical fridge is a purpose-built vaccine fridge designed and built in Australia in line with the National Vaccine Storage guidelines "Strive for 5".

Medisafe Plus vaccine fridges are feature packed, resulting in a vaccine storage solution that gives you peace of mind – knowing your vaccines are stored in the safest possible manner.

The Medisafe Plus 80 is QCPP accredited.

The cloud monitoring acts as an automated data logger which allows you to view and print your temperature data in real time, from any location, on any web enabled smartphone, tablet, or computer.

All units also the ability to have a BMS system plugged into the back of the unit for hospital remote monitoring.

Features

- 24 months Parts and Labour warranty
- Dedicated 24-hour technical support
- QCPP cold chain accredited
- Min/Max Controller
- Temperature Controller
- Dual temperature probes for ambient and simulated vaccine temperatures
- AQ Box Cloud Monitoring – with SMS and email alerts of power or temperature failure
- Audio/Visual alarms
- Lock and Key
- Adjustable shelves
- Low Energy Consumption

Optional Extras

- UPS Unit with Back up battery
- Preventative Maintenance program with extended warranty

DIMENSIONS (HxWxD)	830 x 500 x 560
GROSS WEIGHT	53 Kg
GROSS CAPACITY	80 Lt.
TEMPERATURE RANGE	2° to 8°C
NUMBER OF SHELVES	3
POWER CONSUMPTION PER 24Hr	3.06 kWh
REFRIGERANT	R600a

Local distributor:



SURGICAL AND MEDICAL SUPPLIES PTY LTD
 16 Kensington Road Rose Park SA 5067
 Telephone: (08) 8332 1666 Fax: (08) 8332 5132
 Email: sales@surmed.com.au
 Website: www.surmed.com.au