



Supplied in SA by:



**SURGICAL AND MEDICAL SUPPLIES PTY LTD**  
16 Kensington Road Rose Park SA 5067  
Telephone: (08) 8332 1666 Fax: (08) 8332 5132  
Email: sales@surmed.com.au  
Website: www.surmed.com.au



PHARMA 3000



# PHARMA 3000

FLOOR STANDING  
VACCINE REFRIGERATOR

SPECIFICATIONS



## TEMPERATURE RANGE

+3°C to +7°C  
@ 10°C to 32°C room temperature  
& 60% humidity



## CAPACITY

Effective capacity: 350 Litres



## DIMENSIONS

Width: 595mm  
Depth: 640mm  
Height: 1900mm

FEATURES



## FEATURES

- Purpose built for vaccine storage
- Exceptional performance on all shelves
- Full use of internal space (except 4cm below the fan)
- QCPP approved and exceeds requirements
- No frozen back wall and no cold spots
- Simple to use and no setup required
- Very easy to use hi/lo/reset temperature display
- Can be built in – only a 3cm top gap is required
- Long lasting LED lighting (glass door model)
- 9 adjustable shelves plus raised base shelf
- Designed and manufactured in Denmark
- Lockable and self-closing, glass or solid door
- Available with right or left hinged door
- Energy consumption – only 1.1 kWh/24h
- Extremely quiet – only 36dB
- Audible and visual alarm
- Voltage free external alarm contacts for BMS
- Data logger access valve
- Environmentally responsible refrigerant
- Built in door alarm