

Medisafe Biomedical 98ULT

Medisafe Biomedical 98ULT

Is an undercounter ultra low temperature unit in the Medisafe Biomedical range.

Medisafe Biomedical 98ULT medical freezer is a purpose-built medical freezer designed for the holding of samples and vaccines that require ultra low temp holding

Medisafe Biomedical fridges/freezers are feature packed, resulting in a storage solution that gives you peace of mind – knowing your vaccines/drugs/samples are stored in the safest possible manner.

The Medisafe Biomedical 98ULT is TGA approved.

The optional cloud monitoring acts as an automated data logger which allows you to view and print your temperature data in real time, from any location, on any web enabled smartphone, tablet, or computer.

All units have BMS system plugged into the back of the unit for hospital remote monitoring.



Features

- 24 months Parts and Labour warranty
- Dedicated 24-hour technical support
- Temperature Controller
- Audio/Visual/Probe/ Open Door/Power Failure alarms
- Lock and Key
- S/S shelves
- Battery Back Up for Controller
- Portholes
- LED lighting

Optional Extras

- UPS Unit with Back up battery
- Preventative Maintenance program with extended warranty
- AQ Box Cloud Monitoring– with SMS and email alerts of power and temperature breaches

DIMENSIONS(HxWxD)	825 x 595x 645
GROSS WEIGHT	61 Kg
GROSS CAPACITY	92 Lt.
TEMPERATURE RANGE	-60° to -86°C
NUMBER OF SHELVES	3
POWER CONSUMPTION PER 24Hr	6.7 kWh
REFRIGERANT	Nature R

Local distributor:

Medisafe
biomedical



SURGICAL AND MEDICAL SUPPLIES PTY LTD
16 Kensington Road Rose Park SA 5067
Telephone: (08) 8332 1666 Fax: (08) 8332 5132
Email: sales@surmed.com.au
Website: www.surmed.com.au