



**SURGICAL AND MEDICAL SUPPLIES PTY LTD**

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

**Email: [sales@surmed.com.au](mailto:sales@surmed.com.au)**

**Website: [www.surmed.com.au](http://www.surmed.com.au)**

**INDEX**

	Page
Blanket Warmers	2
Fluid Warmers	4






SURGICAL & MEDICAL SUPPLIES PTY LTD



# Blanket Warmers

## Features:

- Fan forced system
- TGA Approved
- Overheating alarm
- Double glazed doors
- Stainless int./ext. with chrome plated shelves
- Micro-processor control
- Temperature lock out
- Digital Display with Built In Alarms
- Hi / Lo Read out



60 °C



80 L  
257 L  
365 L



500w x 575d x 600h  
850w x 575d x 900h  
750w x 700d x 1800h



## Blanket Warmers



Model	BW5	BW10	BW21
Capacity	80 L	257 L	365 L
External Dimensions	500w x 575d x 600h	850w x 575d x 900h	750w x 700d x 1800h
Internal Dimensions	390w x 480d x 430h	740w x 485d x 590h	630w x 610d x 1590h
Weight	75. Kg	125 Kg	175 Kg
Voltage	240 V	240 V	240 V
Warranty	2 Years general, 3 Years compressor*		
Shelves size	390 x 480	740 x 485	630 x 610
Number of Shelves	2	3	5
Digital Display	✓	✓	✓
Castors	X	✓	✓
Interior/Exterior Cabinet Finish	Stainless Steel	Stainless Steel	Stainless Steel
High/Low Temperature Alarm	✓	✓	✓
Glass Door	✓	✓	✓

### Combination Blanket and Fluid Warmers

Model	Configuration
BW5/FW5	units are stacked
BW10/FW10	units are stacked
BW21/FW21	side-by-side units

\*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.




SURGICAL & MEDICAL SUPPLIES PTY LTD



# Fluid Warmers

## Features:

- Fan forced system
- TGA Approved
- Overheating alarm
- Double glazed doors
- Stainless int./ext. with chrome plated shelves
- Micro-processor control
- Temperature lock out
- Digital Display with Built In Alarms
- Hi / Lo Read out



35 °C



80 L

257 L

365 L



500w x 575d x 600h

850w x 575d x 900h

750w x 700d x 1800h



## Fluid Warmers



Model	FW5	FW10	FW21
Capacity	80 L	257 L	365 L
External Dimensions	500w x 575d x 600h	850w x 575d x 900h	750w x 700d x 1800h
Internal Dimensions	390w x 480d x 430h	740w x 485d x 590h	630w x 610d x 1590h
Weight	75. Kg	125 Kg	175 Kg
Voltage	240 V	240 V	240 V
Warranty	2 Years general, 3 Years compressor*		
Shelves size	390 x 480	740 x 485	630 x 610
Number of Shelves	2	3	5
Digital Display	✓	✓	✓
Castors	X	✓	✓
Interior/Exterior Cabinet Finish	Stainless Steel	Stainless Steel	Stainless Steel
High/Low Temperature Alarm	✓	✓	✓
Glass Door	✓	✓	✓

### Combination Blanket and Fluid Warmers

Model	Configuration
BW5/FW5	units are stacked
BW10/FW10	units are stacked
BW21/FW21	side-by-side units

\*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.