



SURGICAL & MEDICAL SUPPLIES PTY LTD

## Stainless Steel Trolleys



Instrument Trolleys	2
Dressing Trolleys	3
Powder Coated Trolleys	6
Resuscitation Trolleys	7
Other Trolleys	8
Trolley Accessories	8

For enquiries, contact: Chris Lee, Product Specialist  
m. 0476 777 767 e. cleee@surmed.com.au



**SURGICAL AND MEDICAL SUPPLIES PTY LTD**  
16 Kensington Road Rose Park SA 5067  
Telephone: (08) 8332 1666 Fax: (08) 8332 5132  
Email: [sales@surmed.com.au](mailto:sales@surmed.com.au)  
Website: [www.surmed.com.au](http://www.surmed.com.au)

# Instrument Trolleys

## Instrument Trolley - No Rails



- 2 Shelves, no rails
- Full under bracing support for all shelves providing maximum stability & strength
- Double safety folds to all shelf edges
- 2 locking castors
- 304 Grade stainless steel construction - Satin finish

Code	Width	Depth	Height*	Castors
GD0870	490mm	490mm	900mm	100mm
SSTNR50	500mm	500mm	900mm	95mm
GD0880	600mm	490mm	900mm	100mm
SSTNR60	600mm	500mm	900mm	95mm
GD1040	750mm	490mm	900mm	100mm
SSTNR80	800mm	500mm	900mm	95mm
GD0890	900mm	490mm	900mm	100mm
GD0920	1000mm	490mm	900mm	125mm
GD0930	1000mm	600mm	900mm	125mm
GD0940	1200mm	490mm	900mm	125mm
GD1050	1200mm	600mm	900mm	125mm

\*Height is overall including castors

## Instrument Trolley - With Rails



- 2 Shelves, each with rails

Code	Width	Depth	Height*	Castors
GD0900	490mm	490mm	970mm	100mm
SST00	500mm	500mm	900mm	95mm
GD1060	600mm	490mm	970mm	100mm
1450.6	600mm	500mm	900mm	100mm
GD0910	750mm	490mm	900mm	100mm
SSTD00	800mm	500mm	900mm	95mm
GD1070	900mm	490mm	900mm	100mm
GD1080	1000mm	490mm	970mm	125mm
GD1090	1000mm	600mm	970mm	125mm
1450.11	1100mm	500mm	970mm	100mm
GD0860	1200mm	490mm	970mm	125mm
GD1100	1200mm	600mm	970mm	125mm

\*Height is overall including castors. Height to top shelf is 900mm.

# Dressing Trolleys

## Dressing Trolley - 1 Drawer



- 1 x 125mm Drawer with full extension runners
- Full under bracing support for all shelves providing maximum stability & strength
- Double safety folds to all shelf edges
- Work surface with 3 rails, Lower shelf with 3 rails
- 2 locking castors
- 304 Grade stainless steel construction

Code	Width	Depth	Height*	Castors
GD0840	490mm	490mm	970mm	100mm
SST01	500mm	500mm	970mm	95mm
SST01LK†	500mm	500mm	970mm	95mm
GD0950	600mm	490mm	970mm	100mm
1451.6	600mm	500mm	970mm	100mm
GD0850	750mm	490mm	970mm	100mm
SSTD01	800mm	500mm	970mm	95mm
GD0960	900mm	490mm	970mm	100mm

\*Height is overall including castors. Height to top shelf is 900mm. † Locking Drawer

## Dressing Trolley - 2 Drawers



- 2 x 125mm Drawer with full extension runners

Code	Width	Depth	Height*	Castors
GD0970	490mm	490mm	970mm	100mm
SST02	500mm	500mm	970mm	95mm
SST02LK†	500mm	500mm	970mm	95mm
GD0370	600mm	490mm	970mm	100mm
1453.6	600mm	500mm	970mm	100mm
GD0790	750mm	490mm	970mm	100mm
SSTD02	800mm	500mm	970mm	95mm
GD0980	900mm	490mm	970mm	100mm

\*Height is overall including castors. Height to top shelf is 900mm. † Locking Top Drawer

# Dressing Trolleys

## Dressing Trolley - 3 Drawer



- 3 x 125mm Drawer with full extension runners
- Full under bracing support for all shelves providing maximum stability & strength
- Double safety folds to all shelf edges
- Work surface with 3 rails, Lower shelf with 3 rails
- 2 locking castors
- 304 Grade stainless steel construction

Code	Width	Depth	Height*	Castors
GD0990	490mm	490mm	970mm	100mm
SST03	500mm	500mm	970mm	95mm
SST03LK†	500mm	500mm	970mm	95mm
GD1000	600mm	490mm	970mm	100mm
1456	600mm	500mm	970mm	100mm
GD1010	750mm	490mm	970mm	100mm
1456.8	800mm	500mm	970mm	100mm
GD0800	900mm	490mm	970mm	100mm

\*Height is overall including castors. Height to top shelf is 900mm. † Locking Top Drawer

## Dressing Trolley - 4 Drawers



- 4 x 125mm Drawer with full extension runners

Code	Width	Depth	Height*	Castors
GD0810	490mm	490mm	970mm	100mm
GD0810CL‡	490mm	490mm	970mm	100mm
SST04	500mm	500mm	970mm	95mm
SST04LK†	500mm	500mm	970mm	95mm
GD0820	600mm	490mm	970mm	100mm
GD0830	750mm	490mm	970mm	100mm
GD1030	900mm	490mm	970mm	100mm

\*Height is overall including castors † Locking Top Drawer ‡ Central Locking, all Drawers

# Dressing Trolleys

## Dressing Trolley - 4 Drawer



- 2 x 125mm, 2x 250mm Drawers with full extension runners
- Work surface with 3 rails
- 2 locking castors
- 304 Grade stainless steel construction

Code	Width	Depth	Height*	Castors
SSAT04	500mm	500mm	1070mm	95mm
GD0210	600mm	490mm	970mm	125mm

\*Height is overall including castors. Height to top shelf is 900mm.

## Dressing Trolley - 4 Drawer



- 3 x 125mm, 1x 250mm Drawers with full extension runners

Code	Width	Depth	Height*	Castors
GD0230	600mm	490mm	970mm	125mm
1458	600mm	500mm	970mm	100mm

\*Height is overall including castors. Height to top shelf is 900mm.

## Dressing Trolley - 5 Drawer



- 5 x 125mm Drawers with full extension runners

Code	Width	Depth	Height*	Castors
1455	600mm	500mm	970mm	100mm
AX114†	600mm	500mm	970mm	125mm

\*Height is overall including castors. Height to top shelf is 900mm. † With central locking

# Powder Coated Instrument & Dressing Trolleys

- Economical alternative to Stainless Steel (but not for sterile applications)
- White powder-coated steel
- 2 locking castors
- Fully welded



## Powder Coated Instrument Trolley - With Rails

- 1 x 125mm Drawers with full extension runners

Code	Width	Depth	Height*	Castors
PCT00	500mm	500mm	970mm	95mm

\*Height is overall including castors

## Powder Coated Dressing Trolley - 1 Drawer

- 2 x 125mm Drawers with full extension runners

Code	Width	Depth	Height*	Castors
PCT01	500mm	500mm	970mm	95mm

\*Height is overall including castors



## Powder Coated Dressing Trolley - 2 Drawer

- 2 x 125mm Drawers with full extension runners

Code	Width	Depth	Height*	Castors
PCT02	500mm	500mm	970mm	95mm

\*Height is overall including castors. Height to top shelf is 900mm.

## Powder Coated Dressing Trolley - 3 Drawer

- 3 x 125mm Drawers with full extension runners

Code	Width	Depth	Height*	Castors
PCT03	500mm	500mm	970mm	95mm

\*Height is overall including castors



## Powder Coated Dressing Trolley - 4 Drawer

- 4 x 125mm Drawers with full extension runners

Code	Width	Depth	Height*	Castors
PCT04	500mm	500mm	970mm	95mm

\*Height is overall including castors. Height to top shelf is 900mm.

# Resuscitation Trolleys

## Resuscitation Trolley - 2 Half Drawers



- 2 x 125mm Half width drawers
- Work surface with 3 rails
- 2 locking castors
- 304 Grade stainless steel construction

Code	Width	Depth	Height*	Castors
SSTD22	800mm	500mm	970mm	95mm
AX073-G1	900mm	490mm	970mm	125mm

\*Height is overall including castors. Height to top shelf is 900mm.

## Resuscitation Trolley - 4 Half Drawers



- 4 x 125mm Half width drawers

Code	Width	Depth	Height*	Castors
SSTD04	800mm	500mm	970mm	95mm
GC0600	900mm	490mm	970mm	120mm

\*Height is overall including castors. Height to top shelf is 900mm.

## Resuscitation Trolley - 5 Drawers (4 Half, 1 Full)



- 4 x 125mm Half width drawers, 1x 125mm full drawer

Code	Width	Depth	Height*	Castors
1462	800mm	500mm	970mm	100mm

\*Height is overall including castors. Height to top shelf is 900mm.

## Resuscitation Trolley - 7 Drawers (6 Half, 1 Full)



- 6 x 125mm Half width drawers, 1x 250mm full drawer

Code	Width	Depth	Height*	Castors
SSRT07	800mm	500mm	970mm	95mm

\*Height is overall including castors. Height to top shelf is 900mm.

## Other Trolleys

### Anaesthetic Trolley - 6 Drawers



- 1st Drawer 75mm, 2nd & 3rd Drawers 150mm
- 4th Drawer 75mm, 5th & 6th Drawers 100mm
- 2 x Shelves with rails
- 2 x Locking castors
- 304 Grade stainless steel construction

Code	Width	Depth	Height*	Castors
GD0250	900mm	490mm	970mm	125mm

\*Height is overall including castors. Height to top shelf is 900mm.

### Procedure Cart - 6 Drawers



- 4 x 150mm drawers, 2 x 200mm drawers
- 2 x Locking castors
- 304 Grade stainless steel construction

Code	Width	Depth	Height*	Castors
GC0540	900mm	490mm	970mm	120mm

\*Height is overall including castors. Height to top shelf is 900mm.

## Trolley Accessories

### IV Pole with Trolley Mounting Bracket



- Adjustable height stainless steel IV Pole
- Clamp attaches to square trolley posts or legs

Code
1482

### C-Size Oxygen Cylinder Holder



- Stainless steel cylinder holder hooks over a trolley rail
- 125mmW x 110mmD x 350mmH

Code
1480



# Trolley Accessories

## C-Size Oxygen Cylinder Holder



- Stainless steel cylinder holder
- Clamp attaches to square trolley posts or legs

Code

---

AX123

---

