

Distributed in SA by:



SURGICAL AND MEDICAL SUPPLIES PTY LTD

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

Email: sales@surmed.com.au

Website: www.surmed.com.au

Elite
MEDICAL
#becausewecare

Sterile Processing



For enquiries, contact:
Chris Lee, Product Specialist
Surgical & Medical Supplies
m 0476 777 767
e clee@surmed.com.au

www.elitemedical.com.au



“We believe there is always a better way.”

For over 30 years, Elite Medical has stood for people and products of the highest quality to continually supply great value instruments and consumables for health service providers. We keep it simple giving the best quality, life-saving instruments, combined with responsive time-saving service.

While we are not the largest supplier or owned by a multi-national, we are an Australian family-owned business successfully providing the Australian Medical Industry. We know that we can make a difference, so we strive to achieve our partners' goals and dreams **#becausewecare**.

Why Elite Medical

With our focus on **continuous improvement**, we **look after our customers** seeking long-term partnerships, to conduct every transaction with **honesty and transparency**, **standing accountable** to our stakeholders and ourselves for best results.

Regardless of the size of your business or annual spend, our support will always be second to none, and you will always have a dedicated account manager to support you.

Simply order online, give us a call or join our subscription service. Speak to your account manager for more details.

About Neurein®

Backed by our 30 years of experience in premium customer services, our CSSD brand NEUREIN®, provides innovative and reliable solutions engineered to adapt to existing practices introducing a range of consumables tailored for CSSD professionals.



Customer Enquiries

1800 354 836

sales@elitemedical.com.au



Product Offerings

www.elitemedical.com.au

Sterile Processing

<i>Channel Cleaning Brushes</i>	4
<i>Channel Cleaning Brushes</i>	5
<i>General Cleaning Brushes</i>	6
<i>Endoscopy Cleaning Brushes</i>	7
<i>Side Perforated Trays</i>	8
<i>Standard Wire Mesh Baskets</i>	9
<i>Fine Wire Mesh Baskets</i>	10
<i>Accessories</i>	11
<i>SMS Blue Wrap</i>	12
<i>Face Shields</i>	13
<i>Disposable Gowns</i>	14
<i>The perfect solution for your business</i>	15

Channel Cleaning Brushes



1.7mm dia. x 15mm
20018



2mm dia. x 25mm
20010



2.5mm dia. x 25mm
20011



3mm dia. x 25mm
20013



3mm dia. x 55mm
20050



3.5mm dia. x 25mm
20014



4mm dia. x 25mm
20016

Neurein® Channel Cleaning Brushes

CODE	PRODUCT NAME	DESCRIPTION	QTY
20018	Neurein® Channel Cleaning Brush	1.7mm dia. x 15mm Bristle LGTH, 185mm LGTH, Red Bead	PACK/3
20015	Neurein® Channel Cleaning Brush	2mm dia. x 15mm Bristle LGTH, 220mm LGTH, Red Bead	PACK/3
20010	Neurein® Channel Cleaning Brush	2mm dia. x 25mm Bristle LGTH, 220mm LGTH, White Bead	PACK/3
20019	Neurein® Channel Cleaning Brush	2mm dia. x 25mm Bristle LGTH, 300mm LGTH, Red Bead	PACK/3
20011	Neurein® Channel Cleaning Brush	2.5mm dia. x 25mm Bristle LGTH, 220mm LGTH, Blue Bead	PACK/3
20028	Neurein® Channel Cleaning Brush	2.5mm dia. x 25mm Bristle LGTH, 600mm LGTH, Blue Bead	PACK/3
20013	Neurein® Channel Cleaning Brush	3mm dia. x 25mm Bristle LGTH, 220mm LGTH, Yellow Bead	PACK/3
20020	Neurein® Channel Cleaning Brush	3mm dia. x 25mm Bristle LGTH, 300mm LGTH, White Bead	PACK/3
20050	Neurein® Channel Cleaning Brush	3mm dia. x 55mm Bristle LGTH, 600mm LGTH	PACK/3
20014	Neurein® Channel Cleaning Brush	3.5mm dia. x 25mm Bristle LGTH, 220mm LGTH, Black Bead	PACK/3
20016	Neurein® Channel Cleaning Brush	4mm dia. x 25mm Bristle LGTH, 220mm LGTH, Green Bead	PACK/3
20021	Neurein® Channel Cleaning Brush	4mm dia. x 50mm Bristle LGTH, 300mm LGTH, Blue Bead	PACK/3
20051	Neurein® Channel Cleaning Brush	4mm dia. x 55mm Bristle LGTH, 600mm LGTH	PACK/3

Channel Cleaning Brushes

4mm dia. x 25mm
20054



5mm dia. x 50mm
20022



6mm dia. x 70mm
20023



7mm dia. x 80mm
20025



10.5 dia. x 80mm
20026



25mm dia. x 80mm
20036



Neurein® Channel Cleaning Brushes

CODE	PRODUCT NAME	DESCRIPTION	QTY
20054	Neurein® Channel Cleaning Brush	4mm dia. x 25mm Bristle LGTH, 600mm LGTH, 5x Bristle Intervals	PACK/3
20022	Neurein® Channel Cleaning Brush	5mm dia. x 50mm Bristle LGTH, 300mm LGTH, Green Bead	PACK/3
20024	Neurein® Channel Cleaning Brush	5mm dia. x 50mm Bristle LGTH, 400mm LGTH, Green Bead	PACK/3
20029	Neurein® Channel Cleaning Brush	5mm dia. x 50mm Bristle LGTH, 600mm LGTH, Green Bead	PACK/3
20023	Neurein® Channel Cleaning Brush	6mm dia. x 70mm Bristle LGTH, 300mm LGTH, Black Bead	PACK/3
20025	Neurein® Channel Cleaning Brush	7mm dia. x 80mm Bristle LGTH, 400mm LGTH, Yellow Bead	PACK/3
20030	Neurein® Channel Cleaning Brush	7mm dia. x 80mm Bristle LGTH, 600mm LGTH, Yellow Bead	PACK/3
20026	Neurein® Channel Cleaning Brush	10.5mm dia. x 80mm Bristle LGTH, 500mm LGTH, Red Bead	PACK/3
20031	Neurein® Channel Cleaning Brush	10.5mm dia. x 80mm Bristle LGTH, 600mm LGTH, Red Bead	PACK/3
20048	Neurein® Channel Cleaning Brush	12mm dia. x 90mm Bristle LGTH, 200mm LGTH	PACK/3
20027	Neurein® Channel Cleaning Brush	12.5mm dia. x 80mm Bristle LGTH, 400mm LGTH, Blue Bead	PACK/3
20047	Neurein® Channel Cleaning Brush	14mm dia. x 100mm Bristle LGTH, 300mm LGTH	PACK/3
20035	Neurein® Channel Cleaning Brush	17mm dia. x 70mm Bristle LGTH, 400mm LGTH	PACK/3
20036	Neurein® Channel Cleaning Brush	25mm dia. x 80mm Bristle LGTH, 400mm LGTH	PACK/3

General Cleaning Brushes



Neurein® General Cleaning Brushes

CODE	PRODUCT NAME	DESCRIPTION	QTY
20038	Neurein® Toothbrush Cleaning Brush	180mm LGTH, 3 Row, Stainless Steel Bristles	PACK/3
20043	Neurein® Toothbrush Cleaning Brush	180mm LGTH, Nylon Bristles	PACK/3
20045	Neurein® Toothbrush Cleaning Brush	180mm LGTH, Double Ended, 3 Row/1 Row, Nylon Bristles	PACK/3
20074	Neurein® Toothbrush Cleaning Brush	180mm LGTH, 3 Row, Soft Stainless Steel Bristles	PACK/3
20075	Neurein® Toothbrush Cleaning Brush	180mm LGTH, 1 Row, Abrasive Nylon Bristles	PACK/3
20083	Neurein® Toothbrush Cleaning Brush	180mm LGTH, 3 Row, Short Soft Stainless Steel Bristle	PACK/3

Endoscopy Cleaning Brushes

6mm dia. / 6mm dia. x 2300mm
NR-EBF-1004

11mm dia. / 6mm dia. x 2300mm
NR-EBF-1011

6mm dia. / 6 mm dia. x 2300mm LGTH
11mm dia. / 5mm dia. x 160mm LGTH
NR-EBF-1024

3mm dia. / 3mm dia. x 2300mm LGTH
11mm dia. / 5mm dia. x 160mm LGTH
NR-EBF-1026

11mm dia. / 5mm dia/. x 160mm LGTH
NR-EBF-1031



Neurein® Endoscopy Cleaning Brushes

CODE	PRODUCT NAME	DESCRIPTION	QTY
NR-EBF-1004	Neurein® Double-Ended Channel Cleaning Brush	6mm dia. / 6mm dia. x 2300mm LGTH, Single-use	BOX/50
NR-EBF-1011	Neurein® Double-Ended Combination Valve Channel Brush	11mm dia. / 6mm dia. x 2300mm LGTH, Single-use	BOX/100
NR-EBF-1024	Neurein® Double-Ended Channel Brush and Valve Brush Kit	6mm dia. / 6mm dia. x 2300mm LGTH, 11mm dia. / 5mm dia. x 160mm LGTH , Blue Plastic Handle, Single-use	BOX/100
NR-EBF-1026	Neurein® Double-Ended Channel Brush and Valve Brush Kit	3mm dia. / 3mm dia. x 1200mm LGTH, 11mm dia. / 5mm dia. x 160mm LGTH, Blue Plastic Handle, Single-use	BOX/50
NR-EBF-1031	Neurein® Double-Ended Valve Brush	11mm dia. / 5mm dia. x 160mm LGTH, Blue Plastic Handle, Single-use	BOX/200

Side Perforated Trays

Side Perforated and wire mesh base design trays are helpful to contain your instruments during sterilisation process, storage and transportation.

- ⦿ Approximately 5mm x 5mm base mesh size
- ⦿ Compatible with steam and gas plasma sterilisation
- ⦿ Smooth rounded surfaces with no sharp edges



DJ-4550
240mm x 250mm x 60mm

Neurein® Side Perforated Basket

CODE	PRODUCT NAME	DESCRIPTION	QTY
DJ-4545	Neurein® Side Perforated Basket	120 x 250 x 50mm, Flat Base w/ Drop Handles	EACH
DJ-4546	Neurein® Side Perforated Basket	120 x 250 x 60mm, Flat Base, w/ Drop Handles	EACH
DJ-4549	Neurein® Side Perforated Basket	240 x 250 x 50mm, Flat Base w/ Drop Handles	EACH
DJ-4550	Neurein® Side Perforated Basket	240 x 250 x 60mm, Flat Base, w/ Drop Handles	EACH
DJ-4553	Neurein® Side Perforated Basket	405 x 250 x 50mm, Flat Base w/ Drop Handles	EACH
DJ-4554	Neurein® Side Perforated Basket	405 x 250 x 60mm, Flat Base, w/ Drop Handles	EACH
DJ-4557	Neurein® Side Perforated Basket	480 x 250 x 50mm, Flat Base w/ Drop Handles	EACH
DJ-4558	Neurein® Side Perforated Basket	480 x 250 x 60mm, Flat Base, w/ Drop Handles	EACH
DJ-4561	Neurein® Side Perforated Basket	540 x 250 x 50mm, Flat Base, w/ Drop Handles	EACH
DJ-4562	Neurein® Side Perforated Basket	540 x 250 x 60mm, Flat Base, w/ Drop Handles	EACH
DJ-4561-W	Neurein® Side Perforated Basket	540 x 300 x 50mm, Flat Base, w/ Drop Handles	EACH
DJ-4562-W	Neurein® Side Perforated Basket	540 x 300 x 60mm, Flat Base, w/ Drop Handles	EACH

Standard Wire Mesh Baskets

Mesh baskets design trays are helpful to contain your instruments during sterilisation process, storage and transportation.

- ⦿ Approximately 5mm x 5mm mesh size
- ⦿ Compatible with steam and gas plasma sterilisation
- ⦿ Smooth rounded surfaces with no sharp edges



DJ-4413
240mm x 250mm x 60mm

Neurein® Standard Wire Mesh Basket

CODE	PRODUCT NAME	DESCRIPTION	QTY
DJ-4412	Neurein® Standard Wire Mesh Basket	240 x 250 x 50mm, w/ Drop Handles	BOX/1
DJ-4413	Neurein® Standard Wire Mesh Basket	240 x 250 x 60mm, w/ Drop Handles,	BOX/1
DJ-4432	Neurein® Standard Wire Mesh Basket	480 x 250 x 50mm, w/ Drop Handles	BOX/1
DJ-4433	Neurein® Standard Wire Mesh Basket	480 x 250 x 60mm, w/ Drop Handles	BOX/1

Fine Wire Mesh Baskets

Micromesh design trays are helpful to contain your small instruments to ensure they don't get lost during sterilisation process, storage and transportation.

- ⦿ Approximately 2mm x 2mm mesh size
- ⦿ Compatible with steam and gas plasma sterilisation
- ⦿ Smooth rounded surfaces with no sharp edges



150mm x 110mm x 20mm
DJ-4350


Neurein® Fine Mesh Basket

CODE	PRODUCT NAME	DESCRIPTION	QTY
DJ-4343	Neurein® Fine Mesh Basket	40 x 40 x 20mm, w/ Lid	BOX/1
DJ-4345	Neurein® Fine Mesh Basket	80 x 40 x 20mm, w/ Lid	BOX/1
DJ-4348	Neurein® Fine Mesh Basket	80 x 80 x 40mm, w/ Lid	BOX/1
DJ-4350	Neurein® Fine Mesh Basket	150 x 110 x 20mm, w/ Lid	BOX/1
DJ-4351	Neurein® Fine Mesh Basket	160 x 40 x 20mm, w/ Lid	BOX/1
DJ-4357	Neurein® Fine Mesh Basket	200 x 100 x 35mm, w/ Lid	BOX/1
DJ-4360	Neurein® Fine Mesh Basket	235 x 155 x 40mm, w/ Lid	BOX/1

Accessories



Key Surgical®
Surgical Tray Corner Protector
CS-9113



Silicone Pin Mats
DJ-4928



Clinipak®
Tray Liners Sheets
TLC009



Interlock® Plastic Tags
INO1968




Interlock® Metal Tags
INO1449

SMS Blue Wrap

SMS Blue is designed to offer flexibility and enable users to pack conventionally or with 'two layers' at once. SMS Blue is suitable for inner and outer draping of medium to large trays and packs.

- ⦿ Fully synthetic
- ⦿ Single Use



Neurein® SMS Blue Wrap

CODE	PRODUCT NAME	DESCRIPTION	QTY
MLSMS7575	Neurein® SMS Blue Wrapping	75cm X 75cm LGTH	BOX/125
MLSMS9090	Neurein® SMS Blue Wrapping	90cm X 90cm LGTH	BOX/100
MLSMS100100	Neurein® SMS Blue Wrapping	100cm X 100cm LGTH	BOX/100
MLSMS120120	Neurein® SMS Blue Wrapping	120cm X 120cm LGTH	BOX/50

Face Shields

Neurein® medical face shields offer the wearer anti-fog direct splash protection.

- ⦿ Foam headband
- ⦿ Latex-free
- ⦿ Easily disposable



230mm
NR-BK1400XL

Face Shields

CODE	PRODUCT NAME	DESCRIPTION	QTY
NR-BK1400XL	Neurein® Face Shields	230mm LGTH, XL, Single-use	BOX/100

Disposable Gowns

Neurein® washroom-disposable gowns have convenient thumb loops to stop the sleeves from rising, reducing the risk of skin contact with chemicals and bodily fluids.

- ⦿ Impermeable
- ⦿ Back ties
- ⦿ Thumb loops
- ⦿ Single Use



XL
NR-BKSG

Surgical Gowns

CODE	PRODUCT NAME	DESCRIPTION	QTY
NR-BKSG	Neurein® Barrier Gown w/Thumb Loops	XL, Single-use	BOX/100



The perfect solution for your business

Our extensive and trusted range of medical products will fulfil your requirements while delivering improved patient outcomes, ease of use and superior value. We carefully source from only the most respected manufacturers worldwide.

Surgical Consumables

- Diathermy Forceps & Vessel Sealers
- Diathermy Pencils, Tips & Accessories
- Smoke Evacuation Consumables
- Laparoscopic Consumables
- Endoscopy Consumables

Surgical Instruments

- Scalpels & Handles
- Scissors
- Forceps
- Suction Tubes
- Retractors & Skin Hooks
- Needle holders
- Chisels and Gouges
- Fibre Optic Cables

Instrument Care & Sterilisation

- Cleaning Brushes
- Instrument Protection Caps
- Tray Tags & Labels
- Baskets & Trays
- Wraps & Packaging
- Quality Monitoring & Testing
- Transportation & Disposal Aids

Patient Care & Surgical Safety

- Pressure Protection Aids
- Patient Positioning Pads
- Hand Immobilisers
- Steps Stools & Anti-Fatigue Devices
- Specialised Anaesthetics
- PPE

Elite
MEDICAL
#becausewecare

For any enquiry about the products on this catalogue or any of our extensive range, please contact us:

1800 354 836

sales@elitemedical.com.au

www.elitemedical.com.au

Contact our team to discuss your requirements.
We can visit your facility and work with you to
provide the products you need.
If there are products that your department require
but you do not see them in this catalogue, please
contact us as we have a wealth of knowledge for
sourcing bespoke products and solutions.

Distributed in SA by:



SURGICAL AND MEDICAL SUPPLIES PTY LTD

16 Kensington Road Rose Park SA 5067

Telephone: (08) 8332 1666 Fax: (08) 8332 5132

Email: sales@surmed.com.au

Website: www.surmed.com.au

Elite
MEDICAL
#becausewecare

*For any enquiry about the products on this catalogue or
any of our extensive range, please contact us:*

1800 354 836

sales@elitemedical.com.au

www.elitemedical.com.au